

Datasheet for ABIN5658643 Mammaglobin A ELISA Kit



Overview

Quantity:	96 tests
Target:	Mammaglobin A (SCGB2A2)
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum, Tissue Homogenate	
Analytical Method:	Quantitative	
Detection Method:	Colorimetric	
Specificity:	This assay has high sensitivity and excellent specificity for detection of Secretoglobin Family 2A, Member 2 (SCGB2A2). No significant cross-reactivity or interference between Secretoglobir Family 2A, Member 2 (SCGB2A2) and analogues was observed.	
Sensitivity:	0.061 ng/mL	
Target Details		
Target:	Mammaglobin A (SCGB2A2)	
Alternative Name:	Secretoglobin Family 2A, Member 2 (SCGB2A2 Products)	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5658643 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details		
Background:	Gene Name: Secretoglobin Family 2A, Member 2	
	Gene Aliases: MGB1, MGBA, UGB2, Mammaglobin 1, Mammaglobin A, Uteroglobin 2	
Gene ID:	4250	
UniProt:	Q13296	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, air	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate	
	provided in this kit has been pre-coated with an antibody specific to Secretoglobin Family 2A,	
	Member 2 (SCGB2A2). Standards or samples are then added to the appropriate microtiter plate	
	wells with a biotin-conjugated antibody specific to Secretoglobin Family 2A, Member 2	
	(SCGB2A2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each	
	microplate well and incubated. After TMB substrate solution is added, only those wells that	
	contain Secretoglobin Family 2A, Member 2 (SCGB2A2), biotin-conjugated antibody and	
	enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is	
	terminated by the addition of sulphuric acid solution and the color change is measured	
	spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Secretoglobin	
	Family 2A, Member 2 (SCGB2A2) in the samples is then determined by comparing the O.D. of	
	the samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level	
	Secretoglobin Family 2A, Member 2 (SCGB2A2) were tested 20 times on one plate, respectively	
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level	
	Secretoglobin Family 2A, Member 2 (SCGB2A2) were tested on 3 different plates, 8 replicates ir	
	each plate. CV(%) = SD/meanX100	
	Intra-Assay: CV<10%	
	Inter-Assay: CV<12%	
Restrictions:	For Research Use only	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5658643 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months