

Datasheet for ABIN5659548 TMOD3 ELISA Kit



Overview

Quantity:	96 tests
Target:	TMOD3
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA

Product Details

Sample Type:	Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Tropomodulin 3 (TMOD3). No significant cross-reactivity or interference between Tropomodulin 3 (TMOD3) and analogues was observed.
Sensitivity:	0.108 ng/mL
Target Details	

Target:	TMOD3
Alternative Name:	Tropomodulin 3 (TMOD3 Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5659548 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Background:	Gene Name: Tropomodulin 3	
	Gene Aliases: UTMOD, Ubiquitous Tropomodulin	
Gene ID:	50875	
JniProt:	Q9JHJ0	
Pathways:	Regulation of Actin Filament Polymerization	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, ai	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate	
	provided in this kit has been pre-coated with an antibody specific to Tropomodulin 3 (TMOD3).	
	Standards or samples are then added to the appropriate microtiter plate wells with a biotin-	
	conjugated antibody specific to Tropomodulin 3 (TMOD3). Next, Avidin conjugated to	
	Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB	
	substrate solution is added, only those wells that contain Tropomodulin 3 (TMOD3), biotin-	
	conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme	
	substrate reaction is terminated by the addition of sulphuric acid solution and the color change	
	is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of	
	Tropomodulin 3 (TMOD3) in the samples is then determined by comparing the O.D. of the	
	samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level	
	Tropomodulin 3 (TMOD3) were tested 20 times on one plate, respectively	
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level	
	Tropomodulin 3 (TMOD3) were tested on 3 different plates, 8 replicates in each plate. $CV(\%)$ =	
	SD/meanX100	
	Intra-Assay: CV<10%	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5659548 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months