

Datasheet for ABIN5659552 TNNC1 ELISA Kit



Overview

Quantity:	96 tests
Target:	TNNC1
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	62.5 pg/mL - 4000 pg/mL
Minimum Detection Limit:	62.5 pg/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum, Tissue Homogenate		
Analytical Method:	Quantitative		
Detection Method:	Colorimetric		
Specificity:	This assay has high sensitivity and excellent specificity for detection of Troponin C Type 1, Slow (TNNC1). No significant cross-reactivity or interference between Troponin C Type 1, Slow (TNNC1) and analogues was observed.		
Sensitivity:	27.3 pg/mL		
Target Details			
Target:	TNNC1		
Alternative Name:	Troponin C Type 1, Slow (TNNC1 Products)		

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5659552 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details		
Background:	Gene Name: Troponin C Type 1, Slow	
	Gene Aliases: TNC, TroponinC, Cardiac, Troponin C, slow skeletal and cardiac muscles	
Gene ID:	7134	
UniProt:	P63316	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, ai	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate	
	provided in this kit has been pre-coated with an antibody specific to Troponin C Type 1, Slow	
	(TNNC1). Standards or samples are then added to the appropriate microtiter plate wells with a	
	biotin-conjugated antibody specific to Troponin C Type 1, Slow (TNNC1). Next, Avidin	
	conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated.	
	After TMB substrate solution is added, only those wells that contain Troponin C Type 1, Slow	
	(TNNC1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in	
	color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution	
	and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm.	
	The concentration of Troponin C Type 1, Slow (TNNC1) in the samples is then determined by	
	comparing the O.D. of the samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level	
	Troponin C Type 1, Slow (TNNC1) were tested 20 times on one plate, respectively	
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level	
	Troponin C Type 1, Slow (TNNC1) were tested on 3 different plates, 8 replicates in each plate.	
	CV(%) = SD/meanX100	
	Intra-Assay: CV<10%	
	Inter-Assay: CV<12%	
Restrictions:	For Research Use only	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5659552 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months