

## Datasheet for ABIN5659667 BCMA ELISA Kit



Overview

Quantity:	96 tests
Target:	BCMA (TNFRSF17)
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA

## Product Details

Sample Type:	Plasma, Serum
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17). No significant cross-reactivity or interference between Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17) and analogues was observed.
Sensitivity:	0.101 ng/mL
Target Details	

Target:	BCMA (TNFRSF17)
Alternative Name:	Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17 Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5659667 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details	
Background:	Gene Name: Tumor Necrosis Factor Receptor Superfamily, Member 17
	Gene Aliases: CD269, BCM, BCMA, B-cell maturation protein
Gene ID:	21935
UniProt:	088472
Application Details	
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
	5 % within the expiration date under appropriate storage condition. To minimize extra influence
	on the performance, operation procedures and lab conditions, especially room temperature, air
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	the whole assay is performed by the same operator from the beginning to the end.
Assay Time:	3 h
Plate:	Pre-coated
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
	provided in this kit has been pre-coated with an antibody specific to Tumor Necrosis Factor
	Receptor Superfamily, Member 17 (TNFRSF17). Standards or samples are then added to the
	appropriate microtiter plate wells with a biotin-conjugated antibody specific to Tumor Necrosis
	Factor Receptor Superfamily, Member 17 (TNFRSF17). Next, Avidin conjugated to Horseradish
	Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution
	is added, only those wells that contain Tumor Necrosis Factor Receptor Superfamily, Member
	17 (TNFRSF17), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change
	in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution
	and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm.
	The concentration of Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17) in
	the samples is then determined by comparing the O.D. of the samples to the standard curve.
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
	Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17) were tested 20 times
	on one plate, respectively
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
	Tumor Necrosis Factor Receptor Superfamily, Member 17 (TNFRSF17) were tested on 3
	different plates, 8 replicates in each plate. CV(%) = SD/meanX100
	Intra-Assay: CV<10%
	Inter-Assay: CV<12%

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5659667 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Application Details	
Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months