

## Datasheet for ABIN5659792 UROD ELISA Kit



Overview

Quantity:	96 tests
Target:	UROD
Reactivity:	Rat
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA

## Product Details

Alternative Name:

Sample Type:	Tissue Homogenate	
Analytical Method:	Quantitative	
Detection Method:	Colorimetric	
Specificity:	This assay has high sensitivity and excellent specificity for detection of Uroporphyrinogen Decarboxylase (UROD). No significant cross-reactivity or interference between Uroporphyrinogen Decarboxylase (UROD) and analogues was observed.	
Sensitivity:	0.117 ng/mL	
Target Details		
Target:	UROD	

Orde	r at www.antibodies-online.com   www.antikoerper-online.de   www.anticorps-enligne.fr   www.antibodies-online.cn
	International: +49 (0)241 95 163 153   USA & Canada: +1 877 302 8632   support@antibodies-online.com
	Page 1/3   Product datasheet for ABIN5659792   07/25/2024   Copyright antibodies-online. All rights reserved.

Uroporphyrinogen Decarboxylase (UROD Products)

Target Details		
Background:	Gene Name: Uroporphyrinogen Decarboxylase	
	Gene Aliases: PCT, Uroporphyrinogen III Decarboxylase	
UniProt:	P32362	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, air	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Uroporphyrinogen Decarboxylase (UROD). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Uroporphyrinogen Decarboxylase (UROD). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Uroporphyrinogen Decarboxylase (UROD), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically a a wavelength of 450nm ± 10nm. The concentration of Uroporphyrinogen Decarboxylase (UROD) in the samples is then determined by comparing the O.D. of the samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Uroporphyrinogen Decarboxylase (UROD) were tested 20 times on one plate, respectively Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Uroporphyrinogen Decarboxylase (UROD) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10% Inter-Assay: CV<12%	
Restrictions:	For Research Use only	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5659792 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months