

Datasheet for ABIN5659809

K-RAS ELISA Kit



_					
	W	0	rv	10	W

Quantity:	96 tests
Target:	K-RAS (KRAS)
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS). No significant cross-reactivity or interference between V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS) and analogues was observed.
Sensitivity:	0.125 ng/mL

Target Details

Target:	K-RAS (KRAS)
Alternative Name:	V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS Products)

Target Details

Target Type:	Viral Protein
Background:	Gene Name: V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog Gene Aliases: C-K-RAS, K-RAS2A, K-RAS2B, K-RAS4A, K-RAS4B, KI-RAS, KRAS1, KRAS2, NS3, RASK2
Gene ID:	3845
UniProt:	P01116

Application Details

_				
Co	m	m	Or	↑† •
しんし			CI	н.

The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5 % within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

Assay Time:

3 h

Plate:

Pre-coated

Protocol:

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS) in the samples is then determined by comparing the O.D. of the samples to the standard curve.

Assay Precision:

Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS) were tested 20 times on one plate, respectively

Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog (KRAS) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100

Application Details

	Intra-Assay: CV<10%
	Inter-Assay: CV<12%
Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
	specimen samples should be assayed in duplicate. Once the procedure has been started, all
	steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at
	4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant
	pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months