

Datasheet for ABIN5662590

AccuBand™ 100 bp DNA Ladder II[Go to Product page](#)**1** Image

Overview

Quantity:	500 µL
Detection Range:	100 bp - 1000 bp
Minimum Detection Limit:	100 bp
Application:	Agarose Gel Electrophoresis (AGE), SDS-PAGE (SDS)

Product Details

Purpose:	DNA ladder that range in size from 100 ~ 1,000 bp.
Brand:	AccuBand™
Characteristics:	<p>AccuBand™ 100 bp DNA Ladder II is composed of 10 individual DNA fragments, presenting 1k, 900, 800, 700 600, 500, 400, 300, 200 and 100 bp sharp bands respectively.</p> <p>AccuBand™ 100 bp DNA Ladder II provides a sufficient amount of DNA for clear observation of all DNA bands ranging from 100 bp to 1 kb, either in agarose gel or in polyacrylamide gel electrophoresis.</p> <p>AccuBand™ 100 bp DNA Ladder II is ready-to-use, containing loading buffer with dual color tracking dyes (orange G and Xylene cyanol FF)</p>

Application Details

Application Notes:	50 µg/ 500 µL
Restrictions:	For Research Use only

Handling

Format:	Liquid
---------	--------

Handling

Storage: RT, 4 °C, -20 °C

Storage Comment: RT for 6 months, 4°C for 12 months, -20°C for 36 months

Images




Image 1.