



[Go to Product page](#)

Datasheet for ABIN5662608

ExcelBand™ All Blue Regular Range Protein Marker

1 Image

Overview

Quantity: 2 x 250 µL

Detection Range: 10 kDa - 180 kDa

Minimum Detection Limit: 10 kDa

Application: SDS-PAGE (SDS)

Product Details

Purpose: protein molecular weight marker that range in size from 10 to 180 kDa.

Brand: ExcelBand™

Characteristics: The ExcelBand All Blue Regular Range Protein Marker is a blue protein standard with 10 pre-stained proteins covering a wide range of molecular weights from 10 to 180 kDa in Tris-Glycine buffer (9 to 170 kDa in Bis-Tris (MOPS) buffer and Bis-Tris (M

Components: Pre-stained protein marker

Application Details

Application Notes: Optimal working dilution should be determined by the investigator.

Comment: Approximately 0.1~0.5 mg/ml of each protein in the buffer (20 mM Tris-phosphate (pH 7.5), 2% SDS, 0.2 mM DTT, 3.6 M urea, and 15% (v/v) glycerol).

Restrictions: For Research Use only

Handling


Format: Liquid

Handling

Storage: 4 °C, -20 °C

Storage Comment: 4°C for 3 months, -20°C for 24 months

Images



SDS-PAGE

Image 1.