

Datasheet for ABIN5674742  
**GBA ELISA Kit**



[Go to Product page](#)

## Overview

Quantity: 96 tests

Target: GBA

Reactivity: Mouse

Method Type: Sandwich ELISA

Detection Range: 0.156 ng/mL - 10 ng/mL

Minimum Detection Limit: 0.156 ng/mL

Application: ELISA

## Product Details

Sample Type: Plasma, Serum, Tissue Homogenate

Analytical Method: Quantitative

Detection Method: Colorimetric

Sensitivity: 0064 ng/mL

Components: Pre-coated, ready to use 96-well strip plate  
Standard (freeze dried)  
Standard Diluent  
Detection Reagent A  
Detection Reagent B  
Assay Diluent A  
Assay Diluent B  
TMB Substhumane

## Product Details

---

Stop Solution  
Wash Buffer(30 x concenthumane)  
Plate sealer for 96 wells  
Instruction manual

---

Material not included:

1. Microplate reader with  $450 \pm 10$ nm filter.
2. Precision single or multi-channel pipettes and disposable tips.
3. Eppendorf Tubes for diluting samples.
4. Deionized or distilled water.
5. Absorbent paper for blotting the microtiter plate.
6. Container for Wash Solution.

## Target Details

---

Target: GBA

---

Alternative Name: Glucosidase Beta, Acid (GbA) ([GBA Products](#))

---

Background: Alternative name: GBA1, GCB, GLUC, Imiglucerase, Alglucerase, Beta-glucocerebrosidase, Glucosylceramidase,  $\beta$ -Glucosidase, D-Glucosyl-N-Acylsphingosine Glucohydrolase

---

Gene ID: 14466

---

UniProt: [P17439](#)

---

Pathways: [Cellular Glucan Metabolic Process](#)

## Application Details

---

Application Notes: Optimal working dilution should be determined by the investigator.

---

Assay Time: 1 - 4.5 h

---

Plate: Pre-coated

---

Protocol:

1. Prepare all reagents, samples and standards
2. Add 100 $\mu$ L standard or sample to each well. Incubate 2 hours at 37°C
3. Aspirate and add 100 $\mu$ L prepared Detection Reagent A. Incubate 1 hour at 37°C
4. Aspirate and wash 3 times
5. Add 100 $\mu$ L prepared Detection Reagent B. Incubate 1 hour at 37°C
6. Aspirate and wash 5 times
7. Add 90 $\mu$ L Substrate Solution. Incubate 15-25 minutes at 37°C
8. Add 50 $\mu$ L Stop Solution. Read at 450nm immediately.

## Application Details

---

Restrictions: For Research Use only

## Handling

---

Storage: 4 °C

Expiry Date: 12 months