

# Datasheet for ABIN5690979

# **IL12 ELISA Kit**



#### Overview

Quantity:	96 tests
Target:	IL12
Reactivity:	Sheep
Method Type:	Sandwich ELISA
Application:	ELISA

Product Details		
Purpose:	Ovine (sheep) IL-12 p70 ELISA Kit.	
Sample Type:	Cell Culture Supernatant, Cell Samples, Plasma, Serum, Tissue Lysate	
Analytical Method:	Quantitative	
Detection Method:	Colorimetric	
Specificity:	This ELISA antibody pair recognizes Ovine (sheep) IL-12 p70.	
Characteristics:	<ul> <li>Strip plates and additional reagents allow for use in multiple experiments</li> <li>Quantitative protein detection</li> </ul>	
	Establishes normal range	
	The best products for confirmation of antibody array data	
Components:	Pre-Coated 96-well Strip Microplate	
	Wash Buffer	
	Stop Solution	
	Assay Diluent(s)	
	Lyophilized Standard	
	Biotinylated Detection Antibody	

### **Product Details**

- Streptavidin-Conjugated HRP
- TMB One-Step Substrate

#### Material not included:

- · Distilled or deionized water
- Precision pipettes to deliver 2 µL to 1 µL volumes
- Adjustable 1-25 µL pipettes for reagent preparation
- 100 µL and 1 liter graduated cylinders
- Tubes to prepare standard and sample dilutions
- · Absorbent paper
- Microplate reader capable of measuring absorbance at 450nm
- Log-log graph paper or computer and software for ELISA data analysis

## **Target Details**

Target:	IL12	
Alternative Name:	IL-12 p70 (IL12 Products)	
Pathways:	JAK-STAT Signaling, TLR Signaling, Cellular Response to Molecule of Bacterial Origin, Regulation of Leukocyte Mediated Immunity, Positive Regulation of Immune Effector Process,	
	Activated T Cell Proliferation, Cancer Immune Checkpoints, Inflammasome	

## **Application Details**

Plate:	Pre-coated		
Protocol:	Prepare all reagents, samples and standards as instructed in the manual.		
	2. Add 100 µL of standard or sample to each well.		
	3. Incubate 2.5 h at RT or O/N at 4 °C.		
	4. Add 100 μL of prepared biotin antibody to each well.		
	5. Incubate 1 h at RT.		
	6. Add 100 $\mu$ L of prepared Streptavidin solution to each well.		
	7. Incubate 45 min at RT.		
	8. Add 100 μL of TMB One-Step Substrate Reagent to each well.		
	9. Incubate 30 min at RT.		
	10. Add 50 μL of Stop Solution to each well.		
	11. Read at 450 nm immediately.		
Restrictions:	For Research Use only		
Handling			
Storage:	4 °C		

$\vdash$	land	lına
	iaria	шц

Expiry Date:

6 months