

[Go to Product page](#)

Datasheet for ABIN5706400

1X PBS Buffer with 1% Casein

Overview

Quantity:	1 L
Application:	ELISA, FLISA, Immunoprecipitation (IP), Western Blotting (WB)

Product Details

Purification:	1X PBS buffer with 1% casein was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Sterility:	Aseptic filtered

Application Details

Application Notes:	<p>Application Note: This product is a ready-to-use 1X solution and no further dilution is required.</p> <p>Western Blot Dilution: User Optimized</p> <p>Immunoprecipitation Dilution: User Optimized</p> <p>FLISA Dilution: User Optimized</p> <p>ELISA Dilution: User Optimized</p>
Comment:	<p>Synonyms: Phosphate buffered saline, Phosphate buffered solution, PBS, Casein, Western Blot buffer, immunoblot</p> <p>Background: Phosphate buffered saline with casein is suitable for multiple applications including biological diluent buffer for antibodies or other biologics. It also may be used as a wash buffer for immunological assays including western blot, immunohistochemistry, immunofluorescence microscopy, and ELISA. Other applications may require detergents or other additional components.</p>
Restrictions:	For Research Use only

Handling

Format:	Liquid
Concentration:	1 X
Buffer:	Buffer: See application note. Stabilizer: None
Preservative:	Without preservative
Storage:	RT,4 °C
Expiry Date:	6 months