



Datasheet for ABIN5706942

Annexin V-BIOT



[Go to Product page](#)

1 Image

Overview

Quantity:	100 tests
Conjugate:	Biotin
Application:	Flow Cytometry (FACS)

Product Details

Characteristics:	ApoScreen® Annexin V-BIOT
------------------	---------------------------


Application Details

Application Notes:	<ul style="list-style-type: none">• Applications: FC - Quality tested , IHC-WM - Reported in literature , ICC - Reported in literature• Working Dilutions: Flow Cytometry FITC, BIOT, PE, APC, and CY5 conjugates 10 L/106 cells For flow cytometry, the suggested use of these reagents is in a final volume of 100 L of ApoScreen® Annexin Binding Buffer (SB Cat. No. 10045-01)
--------------------	---

Restrictions:	For Research Use only
---------------	-----------------------

Handling

Buffer:	100 tests in 1.0 mL of ApoScreen® Annexin V Binding Buffer
Storage:	4 °C,-20 °C
Storage Comment:	Store at 2-8°C



Flow Cytometry

Image 1. Day old BALB/c mouse bone marrow cells were stained with Annexin V-CY5.