



[Go to Product page](#)

Datasheet for ABIN5706944

Annexin V-FITC

1 Image

Overview

Quantity:	100 tests
Conjugate:	FITC
Application:	Flow Cytometry (FACS)

Product Details

Characteristics:	ApoScreen® Annexin V-FITC
------------------	---------------------------

Application Details

Application Notes:	<ul style="list-style-type: none">• Applications: FC - Quality tested , ICC• Working Dilutions: Flow Cytometry FITC and BIOT conjugates 10 L/10⁶ cells For flow cytometry, the suggested use of these reagents is in a final volume of 100 L of ApoScreen® Annexin Binding Buffer (SB Cat. No. 10045-01)
--------------------	--

Restrictions:	For Research Use only
---------------	-----------------------

Handling

Buffer:	100 tests in 1.0 mL of ApoScreen® Annexin V Binding Buffer
Storage:	4 °C,-20 °C
Storage Comment:	Store at 2-8°C




Image 1. Human T lymphocyte cell line Jurkat was stained with ApoScreen® Annexin V-FITC and ApoScreen® 7-Aminoactinomycin D (7-AAD) at varying stages of apoptosis. (Green - Annexin V-FITC; 7-AAD - Red)