

Datasheet for ABIN6061907  
**anti-CAPZA2 antibody (AA 2-286) (FITC)**



[Go to Product page](#)

## Overview

|                      |                                            |
|----------------------|--------------------------------------------|
| Quantity:            | 100 µL                                     |
| Target:              | CAPZA2                                     |
| Binding Specificity: | AA 2-286                                   |
| Reactivity:          | Human                                      |
| Host:                | Rabbit                                     |
| Clonality:           | Polyclonal                                 |
| Conjugate:           | This CAPZA2 antibody is conjugated to FITC |
| Application:         | ELISA                                      |

## Product Details

|               |                                                                              |
|---------------|------------------------------------------------------------------------------|
| Immunogen:    | Recombinant human F-actin-capping protein subunit alpha-2 protein (2-286AA). |
| Isotype:      | IgG                                                                          |
| Purification: | Protein G                                                                    |

## Target Details

|                   |                                                             |
|-------------------|-------------------------------------------------------------|
| Target:           | CAPZA2                                                      |
| Alternative Name: | CAPZA2 ( <a href="#">CAPZA2 Products</a> )                  |
| Gene ID:          | 830                                                         |
| UniProt:          | <a href="#">P47755</a>                                      |
| Pathways:         | <a href="#">Regulation of Actin Filament Polymerization</a> |

## Application Details

|                    |                                                                    |
|--------------------|--------------------------------------------------------------------|
| Application Notes: | Optimal working dilution should be determined by the investigator. |
|--------------------|--------------------------------------------------------------------|

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

## Handling

|         |        |
|---------|--------|
| Format: | Liquid |
|---------|--------|

|         |                                                           |
|---------|-----------------------------------------------------------|
| Buffer: | 0.01 M PBS, pH 7.4, 0.03 % Proclin-300 and 50 % Glycerol. |
|---------|-----------------------------------------------------------|

|               |         |
|---------------|---------|
| Preservative: | ProClin |
|---------------|---------|

|                    |                                                                                                                   |
|--------------------|-------------------------------------------------------------------------------------------------------------------|
| Precaution of Use: | This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only. |
|--------------------|-------------------------------------------------------------------------------------------------------------------|

|          |               |
|----------|---------------|
| Storage: | -20 °C,-80 °C |
|----------|---------------|

|                  |                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------|
| Storage Comment: | Aliquot and store at -20 °C or -80 °C. Avoid exposure to light. Avoid repeated freeze/thaw cycles. |
|------------------|----------------------------------------------------------------------------------------------------|