

Datasheet for ABIN6066207  
**anti-STT3A antibody (AA 39-114) (FITC)**



[Go to Product page](#)

## Overview

|                      |                                           |
|----------------------|-------------------------------------------|
| Quantity:            | 100 µL                                    |
| Target:              | STT3A                                     |
| Binding Specificity: | AA 39-114                                 |
| Reactivity:          | Human                                     |
| Host:                | Rabbit                                    |
| Clonality:           | Polyclonal                                |
| Conjugate:           | This STT3A antibody is conjugated to FITC |
| Application:         | ELISA                                     |

## Product Details

|               |                                                                                                                    |
|---------------|--------------------------------------------------------------------------------------------------------------------|
| Immunogen:    | Recombinant human Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit STT3A protein (39-114AA). |
| Isotype:      | IgG                                                                                                                |
| Purification: | Protein G                                                                                                          |

## Target Details

|                   |                                          |
|-------------------|------------------------------------------|
| Target:           | STT3A                                    |
| Alternative Name: | STT3A ( <a href="#">STT3A Products</a> ) |
| Gene ID:          | 3703                                     |
| UniProt:          | <a href="#">P46977</a>                   |

## Application Details

|                    |                                                                    |
|--------------------|--------------------------------------------------------------------|
| Application Notes: | Optimal working dilution should be determined by the investigator. |
|--------------------|--------------------------------------------------------------------|

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

## Handling

|         |        |
|---------|--------|
| Format: | Liquid |
|---------|--------|

|         |                                                           |
|---------|-----------------------------------------------------------|
| Buffer: | 0.01 M PBS, pH 7.4, 0.03 % Proclin-300 and 50 % Glycerol. |
|---------|-----------------------------------------------------------|

|               |         |
|---------------|---------|
| Preservative: | ProClin |
|---------------|---------|

|                    |                                                                                                                   |
|--------------------|-------------------------------------------------------------------------------------------------------------------|
| Precaution of Use: | This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only. |
|--------------------|-------------------------------------------------------------------------------------------------------------------|

|          |               |
|----------|---------------|
| Storage: | -20 °C,-80 °C |
|----------|---------------|

|                  |                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------|
| Storage Comment: | Aliquot and store at -20 °C or -80 °C. Avoid exposure to light. Avoid repeated freeze/thaw cycles. |
|------------------|----------------------------------------------------------------------------------------------------|