

Datasheet for ABIN6085120

## **anti-RENT1/UPF1 antibody (AA 956-1125) (FITC)**



[Go to Product page](#)

### Overview

|                      |                                                |
|----------------------|------------------------------------------------|
| Quantity:            | 100 µL                                         |
| Target:              | RENT1/UPF1 (UPF1)                              |
| Binding Specificity: | AA 956-1125                                    |
| Reactivity:          | Human                                          |
| Host:                | Rabbit                                         |
| Clonality:           | Polyclonal                                     |
| Conjugate:           | This RENT1/UPF1 antibody is conjugated to FITC |
| Application:         | ELISA                                          |

### Product Details

|               |                                                                             |
|---------------|-----------------------------------------------------------------------------|
| Immunogen:    | Recombinant Human Regulator of nonsense transcripts 1 protein (956-1125AA). |
| Isotype:      | IgG                                                                         |
| Purification: | Protein G                                                                   |

### Target Details

|                   |                                                |
|-------------------|------------------------------------------------|
| Target:           | RENT1/UPF1 (UPF1)                              |
| Alternative Name: | UPF1 ( <a href="#">UPF1 Products</a> )         |
| Gene ID:          | 5976                                           |
| UniProt:          | <a href="#">Q92900</a>                         |
| Pathways:         | <a href="#">SARS-CoV-2 Protein Interactome</a> |

## Application Details

|                    |                                                                    |
|--------------------|--------------------------------------------------------------------|
| Application Notes: | Optimal working dilution should be determined by the investigator. |
|--------------------|--------------------------------------------------------------------|

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

## Handling

|         |        |
|---------|--------|
| Format: | Liquid |
|---------|--------|

|         |                                                           |
|---------|-----------------------------------------------------------|
| Buffer: | 0.01 M PBS, pH 7.4, 0.03 % Proclin-300 and 50 % Glycerol. |
|---------|-----------------------------------------------------------|

|               |         |
|---------------|---------|
| Preservative: | ProClin |
|---------------|---------|

|                    |                                                                                                                   |
|--------------------|-------------------------------------------------------------------------------------------------------------------|
| Precaution of Use: | This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only. |
|--------------------|-------------------------------------------------------------------------------------------------------------------|

|          |               |
|----------|---------------|
| Storage: | -20 °C,-80 °C |
|----------|---------------|

|                  |                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------|
| Storage Comment: | Aliquot and store at -20 °C or -80 °C. Avoid exposure to light. Avoid repeated freeze/thaw cycles. |
|------------------|----------------------------------------------------------------------------------------------------|