

Datasheet for ABIN6107947

**anti-ADAMTS15 antibody (AA 801-913) (FITC)**[Go to Product page](#)

## Overview

|                      |                                              |
|----------------------|----------------------------------------------|
| Quantity:            | 100 µL                                       |
| Target:              | ADAMTS15                                     |
| Binding Specificity: | AA 801-913                                   |
| Reactivity:          | Human                                        |
| Host:                | Rabbit                                       |
| Clonality:           | Polyclonal                                   |
| Conjugate:           | This ADAMTS15 antibody is conjugated to FITC |
| Application:         | ELISA                                        |

## Product Details

|               |                                                                                                          |
|---------------|----------------------------------------------------------------------------------------------------------|
| Immunogen:    | Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 15 protein (801-913AA). |
| Isotype:      | IgG                                                                                                      |
| Purification: | Protein G                                                                                                |

## Target Details

|                   |                                                |
|-------------------|------------------------------------------------|
| Target:           | ADAMTS15                                       |
| Alternative Name: | ADAMTS15 ( <a href="#">ADAMTS15 Products</a> ) |
| Gene ID:          | 170689                                         |
| UniProt:          | <a href="#">Q8TE58</a>                         |

## Application Details

|                    |                                                                    |
|--------------------|--------------------------------------------------------------------|
| Application Notes: | Optimal working dilution should be determined by the investigator. |
|--------------------|--------------------------------------------------------------------|

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

## Handling

|         |        |
|---------|--------|
| Format: | Liquid |
|---------|--------|

|         |                                                           |
|---------|-----------------------------------------------------------|
| Buffer: | 0.01 M PBS, pH 7.4, 0.03 % Proclin-300 and 50 % Glycerol. |
|---------|-----------------------------------------------------------|

|               |         |
|---------------|---------|
| Preservative: | ProClin |
|---------------|---------|

|                    |                                                                                                                   |
|--------------------|-------------------------------------------------------------------------------------------------------------------|
| Precaution of Use: | This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only. |
|--------------------|-------------------------------------------------------------------------------------------------------------------|

|          |               |
|----------|---------------|
| Storage: | -20 °C,-80 °C |
|----------|---------------|

|                  |                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------|
| Storage Comment: | Aliquot and store at -20 °C or -80 °C. Avoid exposure to light. Avoid repeated freeze/thaw cycles. |
|------------------|----------------------------------------------------------------------------------------------------|