ANTIBODIES ONLINE

Datasheet for ABIN612687 Coagulation Factor V ELISA Kit

Image



Overview

1

Quantity:	96 tests
Target:	Coagulation Factor V (F5)
Reactivity:	Human
Method Type:	Sandwich ELISA
Minimum Detection Limit:	0.9 ng/mL
Application:	ELISA

Product Details

The AssayMax Human Factor V ELISA kit is designed for detection of human Factor V in
plasma, urine, milk, and cell culture supernatants
AssayMax
Plasma, Cell Culture Supernatant
Quantitative
Colorimetric
Factor V Microplate: A 96-well polystyrene microplate (12 strips of 8 wells) coated with a
monoclonal antibody against Factor V. Sealing Tapes: Each kit contains 3 pre-cut, pressure-
sensitive sealing tapes that can be cut to fit the format of the individual assay. Factor V
Standard: Human Factor V in a buffered protein base (60 ng, lyophilized). Biotinylated Factor V
Antibody (80x): A 80-fold concentrated biotinylated polyclonal antibody against human Factor V
(100µl). EIA Diluent Concentrate (10x): A 10-fold concentrated buffered protein base (30 ml).
Wash Buffer Concentrate (20x): A 20-fold concentrated buffered surfactant (30 ml, 2 bottles).

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/5 | Product datasheet for ABIN612687 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

	Streptavidin-Peroxidase Conjugate (SP Conjugate): A 100-fold concentrate (80µl). Chromogen
	Substrate: A ready-to-use stabilized peroxidase chromogen substrate tetramethylbenzidine (8
	ml). Stop Solution: A 0.5 N hydrochloric acid to stop the chromogen substrate reaction (12 ml).
Material not included:	Microplate reader capable of measuring absorbance at 450 nm. Pipettes (1-20 µL, 20-200 µL, 200-1000µLand multiple channel). Deionized or distilled reagent grade water

Target Details

Target:	Coagulation Factor V (F5)
Alternative Name:	Factor V (F5 Products)
Background:	Factor V (FV) is an essential cofactor of the blood coagulation cascade and circulates in plasma as a large single-chain glycoprotein (330 kDa). The deduced amino acid sequence consists of 2224 amino acids inclusive of a 28-amino acid leader peptide. During coagulation, it is converted to the active cofactor FVa via limited proteolysis by thrombin, and spliced into a heavy chain (110 KDa) and a light chain (73 KDa) held together non-covalently by calcium. In the presence of a calcium ion and the phospholipid on cell surfaces, FVa and Fxa form the prothrombinase complex which catalyzes the hydrolysis of prothrombin to thrombin. Thrombin in turn cleaves fibrinogen to fibrin which polymerizes to form a clot. FVa is readily inactivated by anticoagulant activated protein C. FV Leiden, a single amino acid mutation, renders FVa resistant to cleavage by activated protein C. It therefore over-produces thrombin, and leads to excess clotting and hereditary thrombophilia disorder. Severe FV deficiency is associated with mild to severe bleeding diathesis.

Application Details

Sample Volume:	50 µL
Assay Time:	< 4 h
Plate:	Pre-coated
Protocol:	This assay employs a quantitative sandwich enzyme immunoassay technique that measures
	Factor V in less than 4 hours. A polyclonal antibody specific for Factor V has been pre-coated
	onto a microplate. Factor V in standards and samples is sandwiched by the immobilized
	antibody and a biotinylated polyclonal antibody specific for Factor V, which is recognized by a
	streptavidin-peroxidase conjugate. All unbound material is then washed away and a peroxidase
	enzyme substrate is added. The color development is stopped and the intensity of the color is
	measured.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/5 | Product datasheet for ABIN612687 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

Application Details

have formed in the concentrate, mix gently until the crystals have completely dissolved. EIA Diluent Concentrate (10x): Dilute EIA Diluent Concentrate 1:10 with reagent grade water. Store for up to 1 month at 2-8°C. Standard Curve: Reconstitute the 60 ng of Human Factor V Standard with 1 ml of EIA Diluent to generate a stock solution of 60 ng/ml. Allow the standard to sit for 10 minutes with gentle agitation prior to making dilutions. Prepare duplicate or triplicate standard points by serially diluting the Factor V standard solution (60 ng/ml) 1:2 with EIA Diluent to produce 30, 15, 7.5, 3.75, 1.875, and 0.938 ng/ml solutions. EIA Diluent serves as the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining	Reagent Preparation:	Freshly dilute all reagents and bring all reagents to room temperature before use. If crystals
for up to 1 month at 2-8°C. Standard Curve: Reconstitute the 60 ng of Human Factor V Standard with 1 ml of EIA Diluent to generate a stock solution of 60 ng/ml. Allow the standard to sit for 10 minutes with gentle agitation prior to making dilutions. Prepare duplicate or triplicate standard points by serially diluting the Factor V standard solution (60 ng/ml) 1:2 with EIA Diluent to produce 30, 15, 7.5, 3.75, 1.875, and 0.938 ng/ml solutions. EIA Diluent serves as the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		have formed in the concentrate, mix gently until the crystals have completely dissolved. EIA
Standard with 1 ml of EIA Diluent to generate a stock solution of 60 ng/ml. Allow the standard to sit for 10 minutes with gentle agitation prior to making dilutions. Prepare duplicate or triplicate standard points by serially diluting the Factor V standard solution (60 ng/ml) 1:2 with EIA Diluent to produce 30, 15, 7.5, 3.75, 1.875, and 0.938 ng/ml solutions. EIA Diluent serves as the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		Diluent Concentrate (10x): Dilute EIA Diluent Concentrate 1:10 with reagent grade water. Store
to sit for 10 minutes with gentle agitation prior to making dilutions. Prepare duplicate or triplicate standard points by serially diluting the Factor V standard solution (60 ng/ml) 1:2 with EIA Diluent to produce 30, 15, 7.5, 3.75, 1.875, and 0.938 ng/ml solutions. EIA Diluent serves as the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		for up to 1 month at 2-8°C. Standard Curve: Reconstitute the 60 ng of Human Factor V
triplicate standard points by serially diluting the Factor V standard solution (60 ng/ml) 1:2 with EIA Diluent to produce 30, 15, 7.5, 3.75, 1.875, and 0.938 ng/ml solutions. EIA Diluent serves as the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		Standard with 1 ml of EIA Diluent to generate a stock solution of 60 ng/ml. Allow the standard
EIA Diluent to produce 30, 15, 7.5, 3.75, 1.875, and 0.938 ng/ml solutions. EIA Diluent serves as the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		to sit for 10 minutes with gentle agitation prior to making dilutions. Prepare duplicate or
the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		triplicate standard points by serially diluting the Factor V standard solution (60 ng/ml) 1:2 with
Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent 30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		EIA Diluent to produce 30, 15, 7.5, 3.75, 1.875, and 0.938 ng/ml solutions. EIA Diluent serves as
30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		the zero standard (0 ng/ml). Any remaining solution should be frozen at -20°C. Standard Point
P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		Dilution [FV] (ng/ml) P1 1 part Standard (60 ng/ml) 60.000 P2 1 part P1 + 1 part EIA Diluent
Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		30.000 P3 1 part P2 + 1 part EIA Diluent 15.000 P4 1 part P3 + 1 part EIA Diluent 7.500 P5 1 part
briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		P4 + 1 part EIA Diluent 3.750 P6 1 part P5 + 1 part EIA Diluent 1.875 P7 1 part P6 + 1 part EIA
solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		Diluent 0.938 P8 EIA Diluent 0.000 Biotinylated Factor V Antibody (80x): Spin down the antibody
Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		briefly and dilute the desired amount of the antibody 1:80 with EIA Diluent. Any remaining
briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining		solution should be frozen at -20°C. Wash Buffer Concentrate (20x): Dilute Wash Buffer
		Concentrate 1:20 with reagent grade water. SP Conjugate (100x): Spin down the SP Conjugate
solution should be frozen at -20°C		briefly and dilute the desired amount of the conjugate 1:100 with EIA Diluent. Any remaining
		solution should be frozen at -20°C.

Sample Collection:Plasma: Collect plasma using one-tenth volume of 0.1 M sodium citrate as an anticoagulant.
Centrifuge samples at 2000 x g for 10 minutes. Dilute samples 1:800 with EIA Diluent. Store the
remaining samples at -20°C or below. Avoid repeated freeze-thaw cycles. (EDTA or Heparin can
also be used as anticoagulant.) Cell Culture Supernatants: Centrifuge cell culture media at 2000
x g for 10 minutes to remove debris. Collect supernatants and assay. Store the remaining
samples at -20°C or below. Avoid repeated freeze-thaw cycles. Urine: Collect urine using sample
pot. Centrifuge samples at 600 x g for 10 minutes. Urine dilution is suggested at 1:2 in EIA
Diluent. Store samples at -20°C or below for up to 3 months. Avoid repeated freeze-thaw cycles.
Milk: Collect milk using sample tube. Centrifuge samples at 600 x g for 10 minutes. Milk dilution
is suggested at 1:2 in EIA Diluent. Store samples at -20°C or below for up to 3 months. Avoid repeated freeze-thaw cycles.

Assay Procedure: Prepare all reagents, working standards and samples as instructed. Bring all reagents to room temperature before use. The assay is performed at room temperature (20 - 30 °C). Remove excess microplate strips from the plate frame and return them immediately to the foil pouch with desiccant inside. Reseal the pouch securely to minimize exposure to water vapor and store in a vacuum desiccator. Add 50 µL of Standard or sample per well. Cover wells and incubate for two hours. Start the timer after the last sample addition. Wash five times with 200 µL of Wash

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/5 | Product datasheet for ABIN612687 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

	Buffer manually. Invert the plate each time and decant the contents, hit it 4-5 times on
	absorbent paper towel to completely remove the liquid. If using a machine wash six times with
	300 μL of Wash Buffer and then invert the plate, decant the contents, hit it 4-5 times on
	absorbent paper towel to completely remove the liquid. Add 50 μL of Biotinylated Factor V
	Antibody to each well and incubate for one hour. Wash the microplate as described above. Add
	50 μ L of Streptavidin-Peroxidase Conjugate per well and incubate for 30 minutes. Turn on the
	microplate reader and set up the program in advance. Wash the microplate as described above.
	Add 50 μ L of Chromogen Substrate per well and incubate for about 20 minutes or till the
	optimal blue color density develops. Gently tap the plate to ensure thorough mixing and break
	the bubbles in the well with pipette tip. Add 50 μL of Stop Solution to each well. The color will
	change from blue to yellow. Read the absorbance on a microplate reader at a wavelength of
	450 nm immediately. If wavelength correction is available, subtract readings at 570 nm from
	those at 450 nm to correct optical imperfections. Otherwise, read the plate at 450 nm only.
	Please note that some unstable black particles may be generated at high concentration points
	after stopping the reaction for about 10 minutes, which will reduce the readings.
Calculation of Results:	Calculate the mean value of the duplicate or triplicate readings for each standard and sample.
	To generate a Standard Curve, plot 4-parameter graph using the standard concentrations on
	the x-axis and the corresponding mean 450 nm absorbance on the y-axis. T The best-fit line can
	be determined by regression analysis using log-log or four-parameter logistic curve-fit.
	Determine the unknown sample concentration from the Standard Curve and multiply the value
	by the dilution factor. Standard Curve The curve is provided for illustration only. A standard
	curve should be generated each time the assay is performed.
Assay Precision:	Intra-assay and inter-assay coefficients of variation were 4.8 % and 7.4% respectively.
Restrictions:	For Research Use only
Handling	
Handling Advice:	Prepare all reagents (working diluent buffer, wash buffer, standards, biotinylated- antibody, and
	SP conjugate) as instructed, prior to running the assay. Prepare all samples prior to running the
	assay. The dilution factors for the samples are suggested in this protocol. However, the user
	should determine the optimal dilution factor. Spin down the SP conjugate vial and the
	biotinylated-antibody vial before opening and using contents. The kit should not be used
	beyond the expiration date.
Storage:	4 °C/-20 °C
Storage Comment:	Store components of the kit at 2-8°C or -20°C upon arrival up to the expiration date. Store SP

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/5 | Product datasheet for ABIN612687 | 07/26/2024 | Copyright antibodies-online. All rights reserved. Conjugate and Biotinylated Antibody at -20°C Store Microplate, Diluent Concentrate (10x), Wash Buffer, Stop Solution, and Chromogen Substrate at 2-8°C Opened unused microplate wells may be returned to the foil pouch with the desiccant packs. Reseal along zip-seal. May be stored for up to 1 month in a vacuum desiccator. Diluent (1x) may be stored for up to 1 month at 2-8°C. Store Standard at 2-8°C before reconstituting with Diluent and at -20°C after reconstituting with Diluent

Images

