

## Datasheet for ABIN6202330

### FLASH CLIA Kit



[Go to Product page](#)

#### Overview

|                          |                        |
|--------------------------|------------------------|
| Quantity:                | 96 tests               |
| Target:                  | FLASH (CASP8AP2)       |
| Reactivity:              | Rat                    |
| Method Type:             | Sandwich ELISA         |
| Detection Range:         | 0.312 ng/mL - 20 ng/mL |
| Minimum Detection Limit: | 0.312 ng/mL            |
| Application:             | ELISA                  |

#### Product Details

|                    |                                                                                                                                            |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Purpose:           | Rat Caspase 8 Associated Protein 2 (CASP8AP2) CLIA Kit is a Sandwich CLIA Kit for use with Tissue homogenates and other biological fluids. |
| Sample Type:       | Tissue Homogenate                                                                                                                          |
| Analytical Method: | Quantitative                                                                                                                               |
| Detection Method:  | Chemiluminescent                                                                                                                           |
| Sensitivity:       | < 0.136 ng/mL                                                                                                                              |

#### Target Details

|                   |                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------|
| Target:           | FLASH (CASP8AP2)                                                                                  |
| Alternative Name: | Caspase 8 Associated Protein 2 (CASP8AP2) ( <a href="#">CASP8AP2 Products</a> )                   |
| Pathways:         | <a href="#">Chromatin Binding</a> , <a href="#">Positive Regulation of Endopeptidase Activity</a> |

## Application Details

---

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes: | <p>Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.</p> <p>Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.</p> <p>Standard Form: Lyophilized</p> |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

---

|        |            |
|--------|------------|
| Plate: | Pre-coated |
|--------|------------|

---

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

## Handling

---

|          |             |
|----------|-------------|
| Storage: | 4 °C/-20 °C |
|----------|-------------|

---

|                  |                                                                                       |
|------------------|---------------------------------------------------------------------------------------|
| Storage Comment: | Upon receipt, store the kit according to the storage instruction in the kit's manual. |
|------------------|---------------------------------------------------------------------------------------|

---

|              |          |
|--------------|----------|
| Expiry Date: | 6 months |
|--------------|----------|