

Datasheet for ABIN6204961

**HIV-1 Tat Interactive Protein 2, 30kDa (HTATIP2) ELISA Kit**[Go to Product page](#)

## Overview

Quantity:	96 tests
Target:	HIV-1 Tat Interactive Protein 2, 30kDa (HTATIP2)
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.313 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.313 ng/mL
Application:	ELISA

## Product Details

Purpose:	Human TIP30 ELISA Kit is an ELISA Kit against TIP30.
Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Sensitivity:	< 0.188 ng/mL

## Target Details

Target:	HIV-1 Tat Interactive Protein 2, 30kDa (HTATIP2)
Alternative Name:	TIP30 ( <a href="#">HTATIP2 Products</a> )
Target Type:	Viral Protein
UniProt:	<a href="#">Q9BUP3</a>

## Application Details

---

Application Notes:	<p>Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.</p> <p>Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.</p> <p>Standard Form: Lyophilized</p>
--------------------	---

---

Plate:	Pre-coated
--------	------------

---

Restrictions:	For Research Use only
---------------	-----------------------

---

## Handling

---

Storage:	4 °C/-20 °C
----------	-------------

---

Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.
------------------	---

---

Expiry Date:	6 months
--------------	----------

---