

## Datasheet for ABIN6205196 IL21 Receptor CLIA Kit



[Go to Product page](#)

### Overview

|                          |                        |
|--------------------------|------------------------|
| Quantity:                | 96 tests               |
| Target:                  | IL21 Receptor (IL21R)  |
| Reactivity:              | Human                  |
| Method Type:             | Sandwich ELISA         |
| Detection Range:         | 0.156 ng/mL - 10 ng/mL |
| Minimum Detection Limit: | 0.156 ng/mL            |
| Application:             | ELISA                  |

### Product Details

|                    |  |
|--------------------|--|
| Purpose:           | Human Interleukin 21 Receptor (IL21R) CLIA Kit is a Sandwich CLIA Kit for use with Tissue homogenates and other biological fluids. |
| Sample Type:       | Tissue Homogenate  |
| Analytical Method: | Quantitative   |
| Detection Method:  | Chemiluminescent   |
| Sensitivity:       | < 0.063 ng/mL  |

### Target Details

|                   |  |
|-------------------|--|
| Target:           | IL21 Receptor (IL21R)                    |
| Alternative Name: | IL21R ( <a href="#">IL21R Products</a> ) |
| UniProt:          | <a href="#">Q9HBE5</a>                   |

## Target Details

---

Pathways: [JAK-STAT Signaling](#)

## Application Details

---

Application Notes: Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Standard Form: Lyophilized

---

Plate: Pre-coated

---

Restrictions: For Research Use only

## Handling

---

Storage: 4 °C/-20 °C

---

Storage Comment: Upon receipt, store the kit according to the storage instruction in the kit's manual.

---

Expiry Date: 6 months