

Datasheet for ABIN6206685

PIM1 ELISA Kit



Overview

- Overview	
Quantity:	96 tests
Target:	PIM1
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA
Product Details	
Purpose:	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 ELISA Kit is an ELISA kit against
	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 ELISA Kit is an ELISA kit against Mouse Proto-oncogene serine/threonine-protein kinase pim-1 (Pim1).
Purpose:	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 (Pim1).
Purpose: Sample Type:	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 (Pim1). Plasma, Serum
Purpose: Sample Type: Analytical Method:	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 (Pim1). Plasma, Serum Quantitative
Purpose: Sample Type: Analytical Method:	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 (Pim1). Plasma, Serum Quantitative
Purpose: Sample Type: Analytical Method: Detection Method:	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 (Pim1). Plasma, Serum Quantitative
Purpose: Sample Type: Analytical Method: Detection Method: Target Details	Mouse Proto-oncogene serine/threonine-protein kinase pim-1 (Pim1). Plasma, Serum Quantitative Colorimetric

Application Details

Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less
	than 5 % within the expiration date under appropriate storage conditions. To minimize
	performance fluctuations, operation procedures and lab conditions should be strictly controlled.
	It is also strongly suggested that the whole assay is performed by the same user throughout.
	Recommended dilutions: Optimal dilutions/concentrations should be determined by the end
	user.
Plate:	Pre-coated
Restrictions:	For Research Use only
Handling	
Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.
Expiry Date:	6 months