.-online.com antibodies

Datasheet for ABIN6212222 Golgin A2 (GOLGA2) ELISA Kit



Overview

| Quantity: | 96 tests |
|--------------------------|------------------------|
| Target: | Golgin A2 (GOLGA2) |
| Reactivity: | Human |
| Method Type: | Sandwich ELISA |
| Detection Range: | 0.313 ng/mL - 20 ng/mL |
| Minimum Detection Limit: | 0.313 ng/mL |
| Application: | ELISA |

Product Details

| Purpose: | Human GOLGA2/GM130 ELISA Kit is an ELISA Kit against GOLGA2/GM130. |
|--------------------|--|
| Sample Type: | Plasma, Serum, Tissue Homogenate |
| Analytical Method: | Quantitative |
| Detection Method: | Colorimetric |
| Sensitivity: | < 0.188 ng/mL |

Target Details

| Target: | Golgin A2 (GOLGA2) |
|-------------------|--------------------------------|
| Alternative Name: | GOLGA2/GM130 (GOLGA2 Products) |
| UniProt: | Q08379 |
| Pathways: | SARS-CoV-2 Protein Interactome |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN6212222 | 09/10/2023 | Copyright antibodies-online. All rights reserved.

| Application Details | |
|---------------------|--|
| Application Notes: | Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. Standard Form: Lyophilized |
| Plate: | Pre-coated |
| Restrictions: | For Research Use only |
| Handling | |
| Storage: | 4 °C/-20 °C |
| Storage Comment: | Upon receipt, store the kit according to the storage instruction in the kit's manual. |
| Expiry Date: | 6 months |