

Datasheet for ABIN6215452

CUBN CLIA Kit



\sim				
()\	/e	r\/		٨
() 1	v C.	ı vı	\Box	ΙV

Quantity: 96 tests Target: CUBN Reactivity: Human Method Type: Sandwich ELISA Detection Range: 0.312 ng/mL - 20 ng/mL Minimum Detection Limit: 0.312 ng/mL Application: ELISA		
Reactivity: Human Method Type: Sandwich ELISA Detection Range: 0.312 ng/mL - 20 ng/mL Minimum Detection Limit: 0.312 ng/mL Application: ELISA	Quantity:	96 tests
Method Type:Sandwich ELISADetection Range:0.312 ng/mL - 20 ng/mLMinimum Detection Limit:0.312 ng/mLApplication:ELISA	Target:	CUBN
Detection Range: 0.312 ng/mL - 20 ng/mL Minimum Detection Limit: 0.312 ng/mL Application: ELISA	Reactivity:	Human
Minimum Detection Limit: 0.312 ng/mL Application: ELISA	Method Type:	Sandwich ELISA
Application: ELISA	Detection Range:	0.312 ng/mL - 20 ng/mL
	Minimum Detection Limit:	0.312 ng/mL
Droduct Dataila	Application:	ELISA
Product Details	Product Details	

Purpose:	Human Cubilin CLIA Kit is a Sandwich CLIA Kit for use with Tissue homogenates, cell lysates
	and other biological fluids.
Sample Type:	Cell Lysate, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Chemiluminescent
Sensitivity:	< 0.112 ng/mL

Target Details

Target:	CUBN	
Alternative Name:	Cubilin (CUBN Products)	
Pathways:	Metabolism of Steroid Hormones and Vitamin D, Lipid Metabolism	

Application Details

Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less	
	than 5 % within the expiration date under appropriate storage conditions. To minimize	
	performance fluctuations, operation procedures and lab conditions should be strictly controlled.	
	It is also strongly suggested that the whole assay is performed by the same user throughout.	
	Recommended dilutions: Optimal dilutions/concentrations should be determined by the end	
	user.	
	Standard Form: Lyophilized	
Plate:	Pre-coated	
Restrictions:	For Research Use only	
Handling		
Storage:	4 °C/-20 °C	
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.	
Expiry Date:	6 months	