

Datasheet for ABIN6216040

Retinoic Acid Induced 12 (RAI12) ELISA Kit[Go to Product page](#)

Overview

Quantity:	96 tests
Target:	Retinoic Acid Induced 12 (RAI12)
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.313 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.313 ng/mL
Application:	ELISA

Product Details

Purpose:	Human Dermal Papilla Derived Protein 6 (DERP6) ELISA Kit is an ELISA kit against Dermal Papilla Derived Protein 6 (DERP6).
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Sensitivity:	< 0.19 ng/mL

Target Details

Target:	Retinoic Acid Induced 12 (RAI12)
Alternative Name:	Dermal Papilla Derived Protein 6 (DERP6) (RAI12 Products)
UniProt:	Q8TE02

Application Details

Application Notes: Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Standard Form: Lyophilized

Plate: Pre-coated

Restrictions: For Research Use only

Handling

Storage: 4 °C/-20 °C

Storage Comment: Upon receipt, store the kit according to the storage instruction in the kit's manual.

Expiry Date: 6 months