

## Datasheet for ABIN6219444

### MAPKAP Kinase 3 ELISA Kit

[Go to Product page](#)

#### Overview

|                          |                            |
|--------------------------|----------------------------|
| Quantity:                | 96 tests                   |
| Target:                  | MAPKAP Kinase 3 (MAPKAPK3) |
| Reactivity:              | Human                      |
| Method Type:             | Sandwich ELISA             |
| Detection Range:         | 1.56 ng/mL - 100 ng/mL     |
| Minimum Detection Limit: | 1.56 ng/mL                 |
| Application:             | ELISA                      |

#### Product Details

|                    |                                                                                                                                                                  |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Purpose:           | Human Mitogen Activated Protein Kinase Activated Protein Kinase 3 ELISA Kit is an ELISA Kit against Mitogen Activated Protein Kinase Activated Protein Kinase 3. |
| Sample Type:       | Plasma, Serum                                                                                                                                                    |
| Analytical Method: | Quantitative                                                                                                                                                     |
| Detection Method:  | Colorimetric                                                                                                                                                     |
| Sensitivity:       | 0.94 ng/mL                                                                                                                                                       |

#### Target Details

|                   |                                                                                                                                                                |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Target:           | MAPKAP Kinase 3 (MAPKAPK3)                                                                                                                                     |
| Alternative Name: | Mitogen Activated Protein Kinase Activated Protein Kinase 3 ( <a href="#">MAPKAPK3 Products</a> )                                                              |
| Pathways:         | <a href="#">MAPK Signaling</a> , <a href="#">Neurotrophin Signaling Pathway</a> , <a href="#">Activation of Innate immune Response</a> , <a href="#">Toll-</a> |

## Target Details

---

[Like Receptors Cascades](#)

## Application Details

---

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes: | <p>Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.</p> <p>Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.</p> <p>Standard Form: Lyophilized</p> |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|        |            |
|--------|------------|
| Plate: | Pre-coated |
|--------|------------|

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

## Handling

---

|          |             |
|----------|-------------|
| Storage: | 4 °C/-20 °C |
|----------|-------------|

|                  |                                                                                       |
|------------------|---------------------------------------------------------------------------------------|
| Storage Comment: | Upon receipt, store the kit according to the storage instruction in the kit's manual. |
|------------------|---------------------------------------------------------------------------------------|

|              |          |
|--------------|----------|
| Expiry Date: | 6 months |
|--------------|----------|