

## Datasheet for ABIN6227724

## Cilp2 CLIA Kit



## Overview

Quantity:	96 tests
Target:	Cilp2
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA
Product Details	
Purpose:	Lluman Cartilaga Intermediata Lavar Protein 2 (CIL D2) CLIA Vitia a Candwigh CLIA Vit for use
ι αιρόσο.	Human Cartilage Intermediate Layer Protein 2 (CILP2) CLIA Kit is a Sandwich CLIA Kit for use
, d. posc.	with Tissue homogenates, cell lysates and other biological fluids.
Sample Type:	
	with Tissue homogenates, cell lysates and other biological fluids.
Sample Type:	with Tissue homogenates, cell lysates and other biological fluids.  Cell Lysate, Tissue Homogenate
Sample Type:  Analytical Method:	with Tissue homogenates, cell lysates and other biological fluids.  Cell Lysate, Tissue Homogenate  Quantitative
Sample Type:  Analytical Method:  Detection Method:  Sensitivity:	with Tissue homogenates, cell lysates and other biological fluids.  Cell Lysate, Tissue Homogenate  Quantitative  Chemiluminescent
Sample Type:  Analytical Method:  Detection Method:	with Tissue homogenates, cell lysates and other biological fluids.  Cell Lysate, Tissue Homogenate  Quantitative  Chemiluminescent
Sample Type:  Analytical Method:  Detection Method:  Sensitivity:	with Tissue homogenates, cell lysates and other biological fluids.  Cell Lysate, Tissue Homogenate  Quantitative  Chemiluminescent

## **Application Details**

Application Notes:	Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less
	than 5 % within the expiration date under appropriate storage conditions. To minimize
	performance fluctuations, operation procedures and lab conditions should be strictly controlled.
	It is also strongly suggested that the whole assay is performed by the same user throughout.
	Recommended dilutions: Optimal dilutions/concentrations should be determined by the end
	user.
	Standard Form: Lyophilized
Plate:	Pre-coated
Restrictions:	For Research Use only
Handling	
Storage:	4 °C/-20 °C
Storage Comment:	Upon receipt, store the kit according to the storage instruction in the kit's manual.
Expiry Date:	6 months