

[Go to Product page](#)

Datasheet for ABIN6230005

## Cholesterol Esterase ELISA Kit

### Overview

|                          |                            |
|--------------------------|----------------------------|
| Quantity:                | 96 tests                   |
| Target:                  | Cholesterol Esterase (CEL) |
| Reactivity:              | Mouse                      |
| Method Type:             | Sandwich ELISA             |
| Detection Range:         | 31.2 pg/mL - 2000 pg/mL    |
| Minimum Detection Limit: | 31.2 pg/mL                 |
| Application:             | ELISA                      |

### Product Details

|                    |  |
|--------------------|--|
| Purpose:           | Mouse Bile salt-activated lipase ELISA Kit is an ELISA kit against Mouse Bile salt-activated lipase (Cel). |
| Sample Type:       | Plasma, Serum  |
| Analytical Method: | Quantitative   |
| Detection Method:  | Colorimetric   |

### Target Details

|                   |   |
|-------------------|---|
| Target:           | Cholesterol Esterase (CEL)                                  |
| Alternative Name: | Bile salt-activated lipase ( <a href="#">CEL Products</a> ) |
| Pathways:         | <a href="#">Lipid Metabolism</a>                            |

## Application Details

---

**Application Notes:**                      **Stability:** The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.

**Recommended dilutions:** Optimal dilutions/concentrations should be determined by the end user.

---

**Plate:**    Pre-coated

---

**Restrictions:**                                      For Research Use only

---

## Handling

---

**Storage:**    4 °C/-20 °C

---

**Storage Comment:**                                      Upon receipt, store the kit according to the storage instruction in the kit's manual.

---

**Expiry Date:**    6 months

---