

## Datasheet for ABIN6231394 **PAQR7 ELISA Kit**



[Go to Product page](#)

### Overview

|                          |                         |
|--------------------------|-------------------------|
| Quantity:                | 96 tests                |
| Target:                  | PAQR7                   |
| Reactivity:              | Pig                     |
| Method Type:             | Sandwich ELISA          |
| Detection Range:         | 62.5 pg/mL - 4000 pg/mL |
| Minimum Detection Limit: | 62.5 pg/mL              |
| Application:             | ELISA                   |

### Product Details

|                    |                                                    |
|--------------------|----------------------------------------------------|
| Purpose:           | Pig PAQR7 ELISA Kit is an ELISA Kit against PAQR7. |
| Sample Type:       | Plasma, Serum                                      |
| Analytical Method: | Quantitative                                       |
| Detection Method:  | Colorimetric                                       |
| Sensitivity:       | 37.5 pg/mL                                         |

### Target Details

|                   |                                          |
|-------------------|------------------------------------------|
| Target:           | PAQR7                                    |
| Alternative Name: | PAQR7 ( <a href="#">PAQR7 Products</a> ) |

## Application Details

---

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes: | Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user. |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

---

|        |            |
|--------|------------|
| Plate: | Pre-coated |
|--------|------------|

---

|               |                       |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

---

## Handling

---

|          |             |
|----------|-------------|
| Storage: | 4 °C/-20 °C |
|----------|-------------|

---

|                  |                                                                                       |
|------------------|---------------------------------------------------------------------------------------|
| Storage Comment: | Upon receipt, store the kit according to the storage instruction in the kit's manual. |
|------------------|---------------------------------------------------------------------------------------|

---