

Datasheet for ABIN6234287  
**vps13d CLIA Kit**



[Go to Product page](#)

## Overview

|                          |                        |
|--------------------------|------------------------|
| Quantity:                | 96 tests               |
| Target:                  | vps13d                 |
| Reactivity:              | Mouse                  |
| Method Type:             | Sandwich ELISA         |
| Detection Range:         | 0.156 ng/mL - 10 ng/mL |
| Minimum Detection Limit: | 0.156 ng/mL            |
| Application:             | ELISA                  |

## Product Details

|                    |  |
|--------------------|--|
| Purpose:           | Mouse Vacuolar Protein Sorting 13D (VPS13D) CLIA Kit is a Sandwich CLIA Kit for use with Tissue homogenates, cell lysates and other biological fluids. |
| Sample Type:       | Cell Lysate, Tissue Homogenate   |
| Analytical Method: | Quantitative   |
| Detection Method:  | Chemiluminescent   |
| Sensitivity:       | < 0.057 ng/mL  |

## Target Details

|                   |   |
|-------------------|---|
| Target:           | vps13d  |
| Alternative Name: | Vacuolar Protein Sorting 13D (VPS13D) ( <a href="#">vps13d Products</a> ) |

## Application Details

---

Application Notes:                      Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout. Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Standard Form: Lyophilized

---

Plate:                                        Pre-coated

---

Restrictions:                              For Research Use only

---

## Handling

---

Storage:                                    4 °C/-20 °C

---

Storage Comment:                        Upon receipt, store the kit according to the storage instruction in the kit's manual.

---

Expiry Date:                              6 months

---