

[Go to Product page](#)

Datasheet for ABIN6234490

PIEZO1 ELISA Kit

Overview

| | |
|--------------|----------------|
| Quantity: | 96 tests |
| Target: | PIEZO1 |
| Reactivity: | Rat |
| Method Type: | Sandwich ELISA |
| Application: | ELISA |

Product Details

| | |
|--------------------|---|
| Purpose: | Rat Piezo-type mechanosensitive ion channel component 1 ELISA Kit is an ELISA kit against Rat Piezo-type mechanosensitive ion channel component 1 (Piezo1). |
| Sample Type: | Plasma, Serum |
| Analytical Method: | Quantitative |
| Detection Method: | Colorimetric |

Target Details

| | |
|-------------------|---|
| Target: | PIEZO1 |
| Alternative Name: | Piezo-type mechanosensitive ion channel component 1 (PIEZO1 Products) |

Application Details

| | |
|--------------------|--|
| Application Notes: | Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5 % within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. |
|--------------------|--|

Application Details

It is also strongly suggested that the whole assay is performed by the same user throughout.
Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

| | |
|--------|------------|
| Plate: | Pre-coated |
|--------|------------|

| | |
|---------------|-----------------------|
| Restrictions: | For Research Use only |
|---------------|-----------------------|

Handling

| | |
|----------|-------------|
| Storage: | 4 °C/-20 °C |
|----------|-------------|

| | |
|------------------|---|
| Storage Comment: | Upon receipt, store the kit according to the storage instruction in the kit's manual. |
|------------------|---|

| | |
|--------------|----------|
| Expiry Date: | 6 months |
|--------------|----------|