

Datasheet for ABIN6384728



Overview

| Quantity: | 96 tests |
|--------------------|---|
| Target: | IL17F |
| Reactivity: | Pig |
| Method Type: | Sandwich ELISA |
| Application: | ELISA |
| Product Details | |
| Purpose: | Porcine Immunoquantitative (PCR-Based) IL-17F ELISA Kit for cell culture supernatants, |
| | plasma, and serum samples. |
| Sample Type: | Cell Culture Supernatant, Plasma, Serum |
| Analytical Method: | Semi-Quantitative |
| Detection Method: | qPCR |
| Characteristics: | IQ-ELISAs employ specific capture antibodies coated on a 96-well PCR plate. Standards and |
| | samples are pipetted into the wells, the target protein in the standards and samples binds to |
| | the immobilized antibody. The wells are washed and the detection affinity reagent is added to |
| | the wells where it binds to any captured antigen. The wells are washed, and primers and PCR |
| | master mix are added to each well. The plate is placed into a real time PCR instrument for |
| | cycling and measurement of DNA amplification. The cycle number where amplification is |
| | detected is proportional to the amount of affinity detection reagent that bound to captured |
| | antigen in each well. |
| Components: | IL-17F Microplate (Item A): 96 well PCR plate coated with anti-Porcine IL-17F |
| | |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN6384728 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

| | • Wash Buffer I Concentrate (20x) (Item B): 25 ml of 20x concentrated solution |
|------------------------|--|
| | Standards (Item C): 2 vials of recombinant Porcine IL-17F |
| | • Assay Diluent A (Item D): 30 ml diluent buffer, 0.09% sodium azide as preservative. For |
| | Standard/Sample (serum/plasma) diluent |
| | Assay Diluent B (Item E): 15 ml of 5x concentrated buffer. For Standard/Sample (cell culture medium/urine) diluent |
| | Detection Affinity Reagent for IL-17F (Item F): 2 vials of a 4x concentrated solution of anti- Porcine IL-17F affinity reagent |
| | IQELISA Detection Reagent (Item G): 1.4ml of a 10x concentrated stock |
| | Primer Solution (Item I): 1.7 ml vial PCR Master Mix (Item J): 1.2 ml vial |
| | |
| | Final Wash Buffer (Item L): 10ml vial of 10x concentrated buffer |
| Material not included: | Real-time PCR instrument, Bio-Rad recommended |
| | Precision pipettes to deliver 2 µL to 1 mL volumes |
| | Adjustable 1-25 mL pipettes for reagent preparation |
| | 100 mL and 1 liter graduated cylinders |
| | Absorbent paper |
| | Distilled or deionized water |
| | Log-log graph paper or computer and software for data analysis |
| | Tubes to prepare standard or sample dilutions |
| | |

- Heating block or water bath capable of $80^{\circ}C$

Target Details

| Target: | IL17F |
|---------------------|--|
| Alternative Name: | IL-17F (IL17F Products) |
| UniProt: | F1RKV0 |
| Pathways: | Cellular Response to Molecule of Bacterial Origin, Positive Regulation of Endopeptidase Activity |
| Application Details | |
| Application Notes: | The Immuno-Quantitative ELISA (IQELISA) kits are an innovative assay platform that combines the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. This results in an assay that is simultaneously familiar and cutting edge and enables the use of only 1/10th the sample volume while also providing 10x more sensitivity than a traditional ELISA. |
| Comment: | The Immuno-Quantitative ELISA (IQELISA™) kits are an innovative assay platform that combines the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. Also |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN6384728 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

| | called immuno-PCR, this detection platform results in an assay that is simultaneously familiar and cutting edge. Compared to traditional ELISA, IQELISA™ enables the use of only 1/10th the |
|------------------|--|
| Sample Volume: | sample volume while also providing 10x more sensitivity. 25 μL |
| Plate: | Pre-coated |
| Protocol: | Prepare all reagents, samples and standards as instructed Add 25 µL standard or sample to each well. Incubate for 2.5 hours at room temperature or overnight at 4°C Add 25 µL detection affinity reagent to each well. Incubate 1 hour at room temperature |
| | 4. Add 25µL of IQELISA Detection Reagent to each well. Incubate 1 hour at 100m temperature |
| | 5. Add 15µL Primer solution 10µL of PCR master mix to each well |
| | 6. Run real-time PCR |
| Restrictions: | For Research Use only |
| Handling | |
| Storage: | 4 °C,-20 °C,-80 °C |
| Storage Comment: | May be stored for up to 6 months at 2° to 8°C from the date of shipment. Standard |
| | (recombinant protein) should be stored at -20°C or -80°C (recommended at -80°C) after |
| | reconstitution. Opened PCR plate or reagents may be stored for up to 1 month at 2° to 8°C. |
| | Note: the kit can be used within one year if the whole kit is stored at -20°C. Avoid repeated |
| | freeze-thaw cycles. |
| Expiry Date: | 6 months |