

Datasheet for ABIN6385246

Kallikrein 11 IQ-ELISA Kit



Overview

Quantity: 96 tests Target: Kallikrein 11 (KLK11) Reactivity: Human Method Type: Sandwich ELISA Application: ELISA Product Details Purpose: Human Immunoquantitative (PCR-Based) Kallikrein 11 ELISA Kit for plasma, and serum samples. Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washe master mix are added to each well. The plate is placed into a real time.	
Reactivity: Human Method Type: Sandwich ELISA Application: ELISA Product Details Purpose: Human Immunoquantitative (PCR-Based) Kallikrein 11 ELISA Kit for plasma, and serum samples. Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
Method Type: Sandwich ELISA Product Details Purpose: Human Immunoquantitative (PCR-Based) Kallikrein 11 ELISA Kit for plasma, and serum samples. Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed.	
Application: ELISA Product Details Purpose: Human Immunoquantitative (PCR-Based) Kallikrein 11 ELISA Kit for plasma, and serum samples. Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed.	
Purpose: Human Immunoquantitative (PCR-Based) Kallikrein 11 ELISA Kit for plasma, and serum samples. Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
Purpose: Human Immunoquantitative (PCR-Based) Kallikrein 11 ELISA Kit for plasma, and serum samples. Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
plasma, and serum samples. Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
Sample Type: Cell Culture Supernatant, Plasma, Serum Analytical Method: Semi-Quantitative Detection Method: qPCR Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	cell culture supernatants,
Analytical Method: Detection Method: QPCR IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
Detection Method: QPCR IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
Characteristics: IQ-ELISAs employ specific capture antibodies coated on a 96-well PC samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
samples are pipetted into the wells, the target protein in the standard the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are washed	
the immobilized antibody. The wells are washed and the detection at the wells where it binds to any captured antigen. The wells are wash	CR plate. Standards and
the wells where it binds to any captured antigen. The wells are wash	ds and samples binds to
	ffinity reagent is added to
master mix are added to each well. The plate is placed into a real time	ed, and primers and PCR
Thater thin are daded to each well. The place to placed into a real till	ne PCR instrument for
cycling and measurement of DNA amplification. The cycle number w	where amplification is
detected is proportional to the amount of affinity detection reagent t	hat bound to captured
antigen in each well.	
Components: • Kallikrein 11 Microplate (Item A): 96 well PCR plate coated with an	iti-Human Kallikrein 11

- Wash Buffer I Concentrate (20x) (Item B): 25 ml of 20x concentrated solution
- · Standards (Item C): 2 vials of recombinant Human Kallikrein 11
- Assay Diluent A (Item D): 30 ml diluent buffer, 0.09% sodium azide as preservative. For Standard/Sample (serum/plasma) diluent
- Assay Diluent B (Item E): 15 ml of 5x concentrated buffer. For Standard/Sample (cell culture medium/urine) diluent
- Detection Affinity Reagent for Kallikrein 11 (Item F): 2 vials of a 4x concentrated solution of anti-Human Kallikrein 11 affinity reagent
- IQELISA Detection Reagent (Item G): 1.4ml of a 10x concentrated stock
- Primer Solution (Item I): 1.7 ml vial
- PCR Master Mix (Item J): 1.2 ml vial
- PCR Preparation buffer (Item K): 1ml vial of 10x concentrated buffer
- · Final Wash Buffer (Item L): 10ml vial of 10x concentrated buffer

Material not included:

- · Real-time PCR instrument, Bio-Rad recommended
- Precision pipettes to deliver 2 µL to 1 mL volumes
- · Adjustable 1-25 mL pipettes for reagent preparation
- · 100 mL and 1 liter graduated cylinders
- · Absorbent paper
- · Distilled or deionized water
- · Log-log graph paper or computer and software for data analysis
- · Tubes to prepare standard or sample dilutions
- Heating block or water bath capable of 80°C

Target Details

Target:	Kallikrein 11 (KLK11)
Alternative Name:	Kallikrein 11 (KLK11 Products)
Gene ID:	11012
UniProt:	Q9UBX7
Pathways:	Complement System

Application Details	
Application Notes:	The Immuno-Quantitative ELISA (IQELISA) kits are an innovative assay platform that combines
	the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. This results in
	an assay that is simultaneously familiar and cutting edge and enables the use of only 1/10th
	the sample volume while also providing 10x more sensitivity than a traditional ELISA.
Comment:	The Immuno-Quantitative ELISA (IQELISA™) kits are an innovative assay platform that

Application Details

	combines the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. Also
	called immuno-PCR, this detection platform results in an assay that is simultaneously familiar
	and cutting edge. Compared to traditional ELISA, IQELISA™ enables the use of only 1/10th the
	sample volume while also providing 10x more sensitivity.
Sample Volume:	25 μL
Plate:	Pre-coated
Protocol:	1. Prepare all reagents, samples and standards as instructed
	2. Add 25 μL standard or sample to each well. Incubate for 2.5 hours at room temperature or overnight at 4 $^{\circ} C$
	3. Add 25 μL detection affinity reagent to each well. Incubate 1 hour at room temperature
	4. Add 25µL of IQELISA Detection Reagent to each well. Incubate 1 hour
	5. Add 15μL Primer solution 10μL of PCR master mix to each well
	6. Run real-time PCR
Restrictions:	For Research Use only
Handling	
Storage:	4 °C,-20 °C,-80 °C
Storage Comment:	May be stored for up to 6 months at 2° to 8°C from the date of shipment. Standard
	(recombinant protein) should be stored at -20°C or -80°C (recommended at -80°C) after
	reconstitution. Opened PCR plate or reagents may be stored for up to 1 month at 2° to 8°C.
	Note: the kit can be used within one year if the whole kit is stored at -20°C. Avoid repeated
	freeze-thaw cycles.
	Heeze-triaw cycles.
Expiry Date:	6 months