## .-online.com antibodies

# Datasheet for ABIN6385373 ANGPTL7 IQ-ELISA Kit



#### Overview

ernatants,
ndards and
es binds to
t is added to
ers and PCR
ment for
cation is
captured

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN6385373 | 09/10/2023 | Copyright antibodies-online. All rights reserved.

#### Product Details

Components:	<ul> <li>ANGPTL7 Microplate (Item A): 96 well PCR plate coated with anti-Human ANGPTL7</li> <li>Wash Buffer I Concentrate (20x) (Item B): 25 ml of 20x concentrated solution</li> <li>Standards (Item C): 2 vials of recombinant Human ANGPTL7</li> <li>Assay Diluent A (Item D): 30 ml diluent buffer, 0.09% sodium azide as preservative. For Standard/Sample (serum/plasma) diluent</li> <li>Assay Diluent B (Item E): 15 ml of 5x concentrated buffer. For Standard/Sample (cell culture medium/urine) diluent</li> <li>Detection Affinity Reagent for ANGPTL7 (Item F): 2 vials of a 4x concentrated solution of anti- Human ANGPTL7 affinity reagent</li> <li>IQELISA Detection Reagent (Item G): 1.4ml of a 10x concentrated stock</li> <li>Primer Solution (Item I): 1.7 ml vial</li> <li>PCR Master Mix (Item J): 1.2 ml vial</li> <li>PCR Preparation buffer (Item K): 1ml vial of 10x concentrated buffer</li> </ul>
Material not included:	<ul> <li>Real-time PCR instrument, Bio-Rad recommended</li> <li>Precision pipettes to deliver 2 µL to 1 mL volumes</li> <li>Adjustable 1-25 mL pipettes for reagent preparation</li> <li>100 mL and 1 liter graduated cylinders</li> <li>Absorbent paper</li> <li>Distilled or deionized water</li> <li>Log-log graph paper or computer and software for data analysis</li> <li>Tubes to prepare standard or sample dilutions</li> <li>Heating block or water bath capable of 80°C</li> </ul>

#### Target Details

Target:	ANGPTL7
Alternative Name:	ANGPTL7 (ANGPTL7 Products)
Gene ID:	10218
UniProt:	043827

### Application Details

Application Notes:	The Immuno-Quantitative ELISA (IQELISA) kits are an innovative assay platform that combines
	the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. This results in
	an assay that is simultaneously familiar and cutting edge and enables the use of only 1/10th
	the sample volume while also providing 10x more sensitivity than a traditional ELISA.
Comment:	The Immuno-Quantitative ELISA (IQELISA™) kits are an innovative assay platform that
	combines the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. Also

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN6385373 | 09/10/2023 | Copyright antibodies-online. All rights reserved.

	called immuno-PCR, this detection platform results in an assay that is simultaneously familiar and cutting edge. Compared to traditional ELISA, IQELISA™ enables the use of only 1/10th the sample volume while also providing 10x more sensitivity.
Sample Volume:	25 µL
Plate:	Pre-coated
Protocol:	<ol> <li>Prepare all reagents, samples and standards as instructed</li> <li>Add 25 μL standard or sample to each well. Incubate for 2.5 hours at room temperature or overnight at 4°C</li> </ol>
	3. Add 25 µL detection affinity reagent to each well. Incubate 1 hour at room temperature
	<ol> <li>Add 25µL of IQELISA Detection Reagent to each well. Incubate 1 hour</li> <li>Add 15µL Primer solution 10µL of PCR master mix to each well</li> </ol>
	6. Run real-time PCR
Restrictions:	For Research Use only
Handling	
Storage:	4 °C,-20 °C,-80 °C
Storage Comment:	May be stored for up to 6 months at 2° to 8°C from the date of shipment. Standard
	(recombinant protein) should be stored at -20°C or -80°C (recommended at -80°C) after
	reconstitution. Opened PCR plate or reagents may be stored for up to 1 month at 2° to 8°C.
	Note: the kit can be used within one year if the whole kit is stored at -20°C. Avoid repeated
	freeze-thaw cycles.
Expiry Date:	6 months