

Datasheet for ABIN6386014

PSMB8 IQ-ELISA Kit



Go to Product page

| _ | | | | |
|-----|-------------------------|-------|--------|-----|
| () | ve | r\/ | | ۸ / |
| | $^{\prime}$ $^{\prime}$ | : I V | \Box | ٧V |

| Quantity: | 96 tests | |
|--------------------|---|--|
| Target: | PSMB8 | |
| Reactivity: | Human | |
| Method Type: | Sandwich ELISA | |
| Application: | ELISA | |
| Product Details | | |
| Purpose: | Human Immunoquantitative (PCR-Based) LMP7 ELISA Kit for cell culture supernatants, plasma, | |
| | and serum samples. | |
| Sample Type: | Cell Culture Supernatant, Plasma, Serum | |
| Analytical Method: | Semi-Quantitative | |
| Detection Method: | qPCR | |
| Characteristics: | IQ-ELISAs employ specific capture antibodies coated on a 96-well PCR plate. Standards and | |
| | samples are pipetted into the wells, the target protein in the standards and samples binds to | |
| | the immobilized antibody. The wells are washed and the detection affinity reagent is added to | |
| | the wells where it binds to any captured antigen. The wells are washed, and primers and PCR | |
| | master mix are added to each well. The plate is placed into a real time PCR instrument for | |
| | cycling and measurement of DNA amplification. The cycle number where amplification is | |
| | detected is proportional to the amount of affinity detection reagent that bound to captured | |
| | antigen in each well. | |
| Components: | LMP7 Microplate (Item A): 96 well PCR plate coated with anti-Human LMP7 | |

- Wash Buffer I Concentrate (20x) (Item B): 25 ml of 20x concentrated solution
- Standards (Item C): 2 vials of recombinant Human LMP7
- Assay Diluent A (Item D): 30 ml diluent buffer, 0.09% sodium azide as preservative. For Standard/Sample (serum/plasma) diluent
- Assay Diluent B (Item E): 15 ml of 5x concentrated buffer. For Standard/Sample (cell culture medium/urine) diluent
- Detection Affinity Reagent for LMP7 (Item F): 2 vials of a 4x concentrated solution of anti-Human LMP7 affinity reagent
- IQELISA Detection Reagent (Item G): 1.4ml of a 10x concentrated stock
- Primer Solution (Item I): 1.7 ml vial
- PCR Master Mix (Item J): 1.2 ml vial
- PCR Preparation buffer (Item K): 1ml vial of 10x concentrated buffer
- Final Wash Buffer (Item L): 10ml vial of 10x concentrated buffer

Material not included:

- · Real-time PCR instrument, Bio-Rad recommended
- Precision pipettes to deliver 2 µL to 1 mL volumes
- · Adjustable 1-25 mL pipettes for reagent preparation
- · 100 mL and 1 liter graduated cylinders
- · Absorbent paper
- · Distilled or deionized water
- · Log-log graph paper or computer and software for data analysis
- · Tubes to prepare standard or sample dilutions
- Heating block or water bath capable of 80°C

Target Details

| Target: | PSMB8 |
|-------------------|---|
| Alternative Name: | LMP7 (PSMB8 Products) |
| Gene ID: | 5696 |
| Pathways: | Mitotic G1-G1/S Phases, DNA Replication, Synthesis of DNA |

| Application Details | |
|---------------------|---|
| Application Notes: | The Immuno-Quantitative ELISA (IQELISA) kits are an innovative assay platform that combines the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. This results in an assay that is simultaneously familiar and cutting edge and enables the use of only 1/10th the sample volume while also providing 10x more sensitivity than a traditional ELISA. |
| Comment: | The Immuno-Quantitative ELISA (IQELISA™) kits are an innovative assay platform that combines the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. Also |

Application Details

| | called immuno-PCR, this detection platform results in an assay that is simultaneously familia | |
|------------------|---|--|
| | and cutting edge. Compared to traditional ELISA, IQELISA™ enables the use of only 1/10th the | |
| | sample volume while also providing 10x more sensitivity. | |
| Sample Volume: | e: 25 µL | |
| Plate: | Pre-coated | |
| Protocol: | 1. Prepare all reagents, samples and standards as instructed | |
| | 2. Add 25 μL standard or sample to each well. Incubate for 2.5 hours at room temperature or overnight at 4 $^{\circ}C$ | |
| | 3. Add 25 μL detection affinity reagent to each well. Incubate 1 hour at room temperature | |
| | 4. Add 25µL of IQELISA Detection Reagent to each well. Incubate 1 hour | |
| | 5. Add 15µL Primer solution 10µL of PCR master mix to each well | |
| | 6. Run real-time PCR | |
| Restrictions: | For Research Use only | |
| Handling | | |
| Storage: | 4 °C,-20 °C,-80 °C | |
| Storage Comment: | May be stored for up to 6 months at 2° to 8°C from the date of shipment. Standard | |
| | (recombinant protein) should be stored at -20°C or -80°C (recommended at -80°C) after | |
| | reconstitution. Opened PCR plate or reagents may be stored for up to 1 month at 2° to 8°C. | |
| | Note: the kit can be used within one year if the whole kit is stored at -20°C. Avoid repeated | |
| | freeze-thaw cycles. | |
| Expiry Date: | 6 months | |