

## Datasheet for ABIN6700951

# **IL-17A/F Protein**

# 2 Images



### Overview

Quantity:	25 μg
Target:	IL-17A/F
Origin:	Rat
Source:	Escherichia coli (E. coli)
Protein Type:	Recombinant
Application:	SDS-PAGE (SDS)

## **Product Details**

Purpose:	Rat Interleukin-17AF Heterodimer Recombinant Protein
Purification:	Interleukin-17AF Heterodimer purity was determined to be greater than 98% as determined by HpLC, analysis by UV-Spectroscopy at 280nm, and by reducing and non-reducing SDS-pAGE.
Purity:	98,00%
Endotoxin Level:	Measured by LAL is typically ≤ 1 EU/μg protein.
Biological Activity Comment:	The activity is determined by the dose-dependent induction of IL-6 production in cultured mouse NIH 3T3 fibroblasts and is typically 20-30 ng/mL.

# Target Details

Target:	IL-17A/F
Alternative Name:	II17a/II17f (IL-17A/F Products)
Background:	Synonyms: IL17 heterodimer, IL17AF heterodimer, CTLA-8 ML-1 dimer, Interleukin 17AF, Interleukin-17AF heterodimer

Background: Interleukin-17AF (IL-17AF) is a member of the IL-17 family of proteins produced by a subset of T cells, called Th17, following stimulation with IL-23. Since IL-17AF is thought to signal through the IL-17RA receptor, its biological function is similar to that of IL-17A in that it induces the production of a variety of chemokines, in addition to airway neutrophilia. In regard to these functions, IL-17AF has less activity than the IL-17A homodimer but, greater activity than the IL-17F homodimer. Human and rat IL-17AF both show activity on mouse cells. Recombinant rat IL-17AF is a non-glycosylated, disulfide-linked heterodimer, containing one IL-17A subunit and one IL-17F subunit. The dimer has a total of 269 amino acids, with total molecular weight of 30.7 kDa.

UniProt:

Q61453

## **Application Details**

Application Notes: Other: User Optimized

Application\_Note: Interleukin-17AF Heterodimer Recombinant Protein has been tested by SDS-PAGE and biological activity and is suitable as a control for polyclonal or monoclonal anti-Interleukin-17AF Heterodimer in immunological assays.

Suggested\_Applications: Cellular Assay

Other\_Performance\_Data:

Restrictions: For Research Use only

### Handling

Comment:


Format:	Lyophilized
Reconstitution:	Reconstitution_Buffer: Restore with deionized water (or equivalent) Reconstitution_Volume: 25 μL (25-250 μL)
Buffer:	Buffer: 0.1 % Trifluoroacetic acid Stabilizer: None
Preservative:	Without preservative
Storage:	4 °C,-20 °C
Storage Comment:	Store vial at 4° C prior to restoration. Dilute only prior to immediate use. Maintain sterility. This product DOES NOT contain preservative. DO NOT VORTEX. We recommend adding a carrier protein such as HSA or BSA to 0.1% (i.e. 1.0 mg/mL). For best results aliquot contents and freeze at -20° C or colder. Avoid cycles of freezing and thawing. Centrifuge vial before each

opening to dislodge contents from the cap and to clarify if contents are not clear after standing at room temperature.

**Expiry Date:** 

6 months


#### **Images**



#### **SDS-PAGE**

**Image 1.** SDS-PAGE of Rat Interleukin-17AF Heterodimer Recombinant Protein SDS-PAGE of Rat Interleukin-17 Animal Free Recombinant Protein. Lane 1: 1  $\mu$ g Rat IL-17 AF in non-reducing conditions . Lane 2: 1  $\mu$ g Rat IL-17 AF in reducing conditions (+). Lane 3: Molecular weight marker. Rat IL-17AF is a heterodimer with a predicted total MW of 30.7 kDa.

#### Rat IL-17AF Induced Production of IL-6 by NIH 3T3 Cells



#### **SDS-PAGE**

Image 2. SDS-PAGE of Rat Interleukin-17AF Heterodimer Recombinant Protein Bioactivity of Rat Interleukin-17 AF Heterodimer Recombinant Protein. Serial dilutions of Rat IL-17 AF (starting at 2.5 ug/mL) were added to NIH 3T3 cells. After 48 hours, production of mouse IL-6 was measured and the linear portion of the curve was us used to calculate the ED50. The ED50 of Rat IL-17 AF is between 15.3-23 ng/mL. There is no typical range for Rat IL-17AF.