.-online.com antibodies

Datasheet for ABIN6730851 Sonic Hedgehog IQ-ELISA Kit



Overview		
Quantity:	96 tests	
Target:	Sonic Hedgehog (SHH)	
Reactivity:	Human	
Method Type:	DNA-Binding ELISA	
Application:	ELISA	
Product Details		
Purpose:	Human Immunoquantitative (PCR-Based) Sonic Hedgehog N-Terminal ELISA Kit for cell culture supernatants, plasma, and serum samples.	
Sample Type:	Cell Culture Supernatant, Plasma, Serum	
Analytical Method:	Quantitative	
Detection Method:	qPCR	
Specificity:	4.2 pg/ml in Assay Diluent C, and8 pg/ml in Assay Diluent B. (Anticipated minimum sensitivity. Assay may detect lower levels of antigen.)	
Sensitivity:	4.2 pg/mL	
Characteristics:	The Immuno Quantitative Enzyme Linked ImumunoSorbent Assay (IQELISA [™]) is an innovative new assay that combines the specificity and ease of use of an ELISA with the sensitivity of real- time PCR. This results in an assay that is simultaneously familiar and cutting edge and enables the use of lower sample volumes while also providing more sensitivity. The Human Sonic Hedgehog N-Terminal IQELISA [™] Kit is a modified ELISA assay with high sensitivity qPCR readout for the quantitative measurement of Human Sonic Hedgehog N-Terminal in serum,	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN6730851 | 09/10/2023 | Copyright antibodies-online. All rights reserved.

	plasma, and cell culture supernatants. This assay employs an antibody specific for Human Sonic Hedgehog N-Terminal coated on a 96-well PCR plate. Standards and samples are pipetted into the wells and Sonic Hedgehog N-Terminal present in a sample is bound to the wells by the immobilized antibody. The wells are washed and a detection affinity molecule is added to the plates. After washing away unbound detection affinity molecule, primers and a PCR master mix are added to the wells and data is collected using qPCR. Ct values obtained from the qPCR are then used to calculate the amount of antigen contained in each sample, where lower Ct values indicate a higher concentration of antigen.
Components:	 Sonic Hedgehog N-Terminal Microplate (Item A): 96 well PCR plate coated with anti-Human Sonic Hedgehog N-Terminal Wash Buffer I Concentrate (20x) (Item B): 25 ml of 20x concentrated solution Standards (Item C): 2 vials of recombinant Human Sonic Hedgehog N-Terminal Assay Diluent A (Item D): 30 ml diluent buffer, 0.09% sodium azide as preservative. For Standard/Sample (serum/plasma) diluent Assay Diluent B (Item E): 15 ml of 5x concentrated buffer. For Standard/Sample (cell culture medium/urine) diluent Detection Affinity Reagent for Sonic Hedgehog N-Terminal (Item F): 2 vials of a 4x concentrated solution of anti-Human Sonic Hedgehog N-Terminal affinity reagent IQELISA Detection Reagent (Item G): 1.4ml of a 10x concentrated stock Primer Solution (Item I): 1.7 ml vial PCR Master Mix (Item J): 1.2 ml vial PCR Preparation buffer (Item K): 1ml vial of 10x concentrated buffer Final Wash Buffer (Item L): 10ml vial of 10x concentrated buffer
Material not included:	 Real-time PCR instrument, Bio-Rad recommended Precision pipettes to deliver 2 µL to 1 mL volumes Adjustable 1-25 mL pipettes for reagent preparation 100 mL and 1 liter graduated cylinders Absorbent paper Distilled or deionized water Log-log graph paper or computer and software for data analysis Tubes to prepare standard or sample dilutions Heating block or water bath capable of 80°C

Target Details

Target:	Sonic Hedgehog (SHH)
Alternative Name:	Sonic Hedgehog (SHH Products)
Gene ID:	6469

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN6730851 | 09/10/2023 | Copyright antibodies-online. All rights reserved.

Target Details	
UniProt:	Q15465
Pathways:	Hedgehog Signaling, Dopaminergic Neurogenesis, Regulation of Muscle Cell Differentiation,
	Tube Formation, Skeletal Muscle Fiber Development
Application Details	
Application Notes:	The Immuno-Quantitative ELISA (IQELISA) kits are an innovative assay platform that combines
	the specificity and ease of use of an ELISA with the sensitivity of real-time PCR. This results in
	an assay that is simultaneously familiar and cutting edge and enables the use of only 1/10th
	the sample volume while also providing 10x more sensitivity than a traditional ELISA.
Comment:	Features:
	- Excellent for follow up/confirmation of antibody array data
	- Pre-coated 96 well plate
	- Strip plate format minimizes waste
	- Complete kit includes all necessary reagents
	- 10x more sensitive than ELISA
	- Reproducible and reliable
	- Optimized for serum, plasma, urine and cell culture medium
	- 1/10th the sample volume of an ELISA
	Research Applications:
	- Detection of quantitative protein levels in biological fluids
	- Validation of antibody array results
	- Validation of biomarker discovery studies
Plate:	Pre-coated
Protocol:	1. Prepare all reagents, samples and standards as instructed
	 Add 25 µL standard or sample to each well. Incubate for 2.5 hours at room temperature or overnight at 4°C
	3. Add 25 μL detection affinity reagent to each well. Incubate 1 hour at room temperature
	4. Add 25µL of IQELISA Detection Reagent to each well. Incubate 1 hour
	5. Add 15 μ L Primer solution 10 μ L of PCR master mix to each well
	6. Run real-time PCR
Restrictions:	For Research Use only

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/4 | Product datasheet for ABIN6730851 | 09/10/2023 | Copyright antibodies-online. All rights reserved.

Handling		
Storage Comment:	May be stored for up to 6 months at 2° to 8°C from the date of shipment. Standard	
	(recombinant protein) should be stored at -20°C or -80°C (recommended at -80°C) after	
	reconstitution. Opened PCR plate or reagents may be stored for up to 1 month at 2°C to 8°C.	
Expiry Date:	6 months	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/4 | Product datasheet for ABIN6730851 | 09/10/2023 | Copyright antibodies-online. All rights reserved.