



[Go to Product page](#)

## Datasheet for ABIN6951350

### CALML3 ELISA Kit

#### Overview

Quantity:	96 tests
Target:	CALML3
Reactivity:	Human
Method Type:	Sandwich ELISA
Application:	ELISA

#### Product Details

Purpose:	Human CALML3 ELISA Kit.
Sample Type:	Cell Culture Supernatant, Cell Samples, Plasma, Serum, Tissue Lysate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Characteristics:	<ul style="list-style-type: none"> <li>• Strip plates and additional reagents allow for use in multiple experiments</li> <li>• Quantitative protein detection</li> <li>• Establishes normal range</li> <li>• The best products for confirmation of antibody array data</li> </ul>
Components:	<ul style="list-style-type: none"> <li>• Pre-Coated 96-well Strip Microplate</li> <li>• Wash Buffer</li> <li>• Stop Solution</li> <li>• Assay Diluent(s)</li> <li>• Lyophilized Standard</li> <li>• Biotinylated Detection Antibody</li> <li>• Streptavidin-Conjugated HRP</li> <li>• TMB One-Step Substrate</li> </ul>

## Product Details

---

Material not included:	<ul style="list-style-type: none"><li>• Distilled or deionized water</li><li>• Precision pipettes to deliver 2 µl to 1 µl volumes</li><li>• Adjustable 1-25 µl pipettes for reagent preparation</li><li>• 100 µl and 1 liter graduated cylinders</li><li>• Tubes to prepare standard and sample dilutions</li><li>• Absorbent paper</li><li>• Microplate reader capable of measuring absorbance at 450nm</li><li>• Log-log graph paper or computer and software for ELISA data analysis</li></ul>
------------------------	---

## Target Details

---

Target:	CALML3
Alternative Name:	CALML3 ( <a href="#">CALML3 Products</a> )
Gene ID:	810
UniProt:	<a href="#">P27482</a>
Pathways:	<a href="#">Phototransduction</a>

## Application Details

---

Application Notes:	Optimal working dilution should be determined by the investigator.
Plate:	Pre-coated
Protocol:	<ol style="list-style-type: none"><li>1. Prepare all reagents, samples and standards as instructed in the manual.</li><li>2. Add 100 µl of standard or sample to each well.</li><li>3. Incubate 2.5 h at RT or O/N at 4°C.</li><li>4. Add 100 µl of prepared biotin antibody to each well.</li><li>5. Incubate 1 h at RT.</li><li>6. Add 100 µl of prepared Streptavidin solution to each well.</li><li>7. Incubate 45 min at RT.</li><li>8. Add 100 µl of TMB One-Step Substrate Reagent to each well.</li><li>9. Incubate 30 min at RT.</li><li>10. Add 50 µl of Stop Solution to each well.</li><li>11. Read at 450 nm immediately.</li></ol>
Restrictions:	For Research Use only

## Handling

---

Expiry Date:	6 months
--------------	----------