

Datasheet for ABIN6981853

**anti-Ras Protein-Specific Guanine Nucleotide-Releasing Factor 2 (RASGRF2) (AA 20-120) antibody (Alexa Fluor 750)**[Go to Product page](#)

## Overview

Quantity:	100 µL
Target:	Ras Protein-Specific Guanine Nucleotide-Releasing Factor 2 (RASGRF2)
Binding Specificity:	AA 20-120
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	Alexa Fluor 750
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

## Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human GRF2
Isotype:	IgG
Cross-Reactivity:	Mouse
Predicted Reactivity:	Human,Rat,Cow,Sheep
Purification:	Purified by Protein A.

## Target Details

Target:	Ras Protein-Specific Guanine Nucleotide-Releasing Factor 2 (RASGRF2)
Alternative Name:	GRF2 ( <a href="#">RASGRF2 Products</a> )

## Target Details

Background:	Synonyms: GRF2, RAS-GRF2, Ras-specific guanine nucleotide-releasing factor 2, Ras guanine nucleotide exchange factor 2, RASGRF2  Background: Functions as a calcium-regulated nucleotide exchange factor activating both Ras and RAC1 through the exchange of bound GDP for GTP. Preferentially activates HRAS in vivo compared to RRAS based on their different types of prenylation. Functions in synaptic plasticity by contributing to the induction of long term potentiation.
Gene ID:	5924
UniProt:	<a href="#">O14827</a>
Pathways:	<a href="#">Neurotrophin Signaling Pathway</a>

## Application Details

Application Notes:	IF(IHC-P) 1:50-200 IF(IHC-F) 1:50-200 IF(ICC) 1:50-200
Restrictions:	For Research Use only

## Handling

Format:	Liquid
Concentration:	1 µg/µL
Buffer:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	-20 °C
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.
Expiry Date:	12 months