antibodies .- online.com





anti-Ras Protein-Specific Guanine Nucleotide-Releasing Factor 2 (RASGRF2) (AA 20-120) antibody (Cy3)



Go to Product page

()	11/0	K\ /	iew	1
	\cup	'I V/I	$\square \vee \vee$	ı

Quantity:	100 μL	
Target:	Ras Protein-Specific Guanine Nucleotide-Releasing Factor 2 (RASGRF2)	
Binding Specificity:	AA 20-120	
Reactivity:	Mouse	
Host:	Rabbit	
Clonality:	Polyclonal	
Conjugate:	СуЗ	
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))	

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human GRF2
Isotype:	IgG
Cross-Reactivity:	Mouse
Predicted Reactivity:	Human,Rat,Cow,Sheep
Purification:	Purified by Protein A.

Target Details

Target:	Ras Protein-Specific Guanine Nucleotide-Releasing Factor 2 (RASGRF2)
Alternative Name:	GRF2 (RASGRF2 Products)

Target Details

rai got botano		
Background:	Synonyms: GRF2, RAS-GRF2, Ras-specific guanine nucleotide-releasing factor 2, Ras guanine	
	nucleotide exchange factor 2, RASGRF2	
	Background: Functions as a calcium-regulated nucleotide exchange factor activating both Ras	
	and RAC1 through the exchange of bound GDP for GTP. Preferentially activates HRAS in vivo	
	compared to RRAS based on their different types of prenylation. Functions in synaptic plasticity	
	by contributing to the induction of long term potentiation.	
Gene ID:	5924	
UniProt:	014827	
Pathways:	Neurotrophin Signaling Pathway	
Application Details		
Application Notes:	IF(IHC-P) 1:50-200	
	IF(IHC-F) 1:50-200	
	IF(ICC) 1:50-200	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Concentration:	1 μg/μL	
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.	
Preservative:	ProClin	
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.	
Storage:	-20 °C	
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.	
Expiry Date:	12 months	