



[Go to Product page](#)

Datasheet for ABIN6984461
anti-NIPA1 antibody (AA 161-260) (Cy5)

Overview

Quantity:	100 µL
Target:	NIPA1
Binding Specificity:	AA 161-260
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This NIPA1 antibody is conjugated to Cy5
Application:	Immunofluorescence (Paraffin-embedded Sections) (IF (p)), Immunofluorescence (Cultured Cells) (IF (cc))

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human NIPA1
Isotype:	IgG
Cross-Reactivity:	Human
Predicted Reactivity:	Mouse,Rat,Dog
Purification:	Purified by Protein A.

Target Details

Target:	NIPA1
Alternative Name:	NIPA1 (NIPA1 Products)

Target Details

Background: Synonyms: FSP 3, FSP3, Magnesium transporter NIPA1, MGC102724, MGC35570, NIPA 1, NIPA1, NIPA1_HUMAN, Non imprinted in Prader Willi/Angelman syndrome region protein 1, Non-imprinted in Prader-Willi/Angelman syndrome region protein 1, Spastic paraplegia 6 (autosomal dominant), Spastic paraplegia 6 protein, SPG 6, SPG6.

Background: This gene encodes a magnesium transporter that associates with early endosomes and the cell surface in a variety of neuronal and epithelial cells. This protein may play a role in nervous system development and maintenance. Multiple transcript variants encoding different isoforms have been found for this gene. Mutations in this gene have been associated with autosomal dominant spastic paraplegia 6. [provided by RefSeq, Nov 2008]

Gene ID: 123606

UniProt: [Q7RTP0](#)

Application Details

Application Notes: IF(IHC-P) 1:50-200
IF(IHC-F) 1:50-200
IF(ICC) 1:50-200

Restrictions: For Research Use only

Handling

Format: Liquid

Concentration: 1 µg/µL

Buffer: Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.

Preservative: ProClin

Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.

Storage: -20 °C

Storage Comment: Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.

Expiry Date: 12 months