



[Go to Product page](#)

Datasheet for ABIN701927

## anti-Tetraspanin 5 antibody (AA 65-160) (FITC)

### Overview

Quantity:	100 µL
Target:	Tetraspanin 5 (TSPAN5)
Binding Specificity:	AA 65-160
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This Tetraspanin 5 antibody is conjugated to FITC
Application:	Western Blotting (WB), Immunofluorescence (Paraffin-embedded Sections) (IF (p)), Immunofluorescence (Cultured Cells) (IF (cc))

### Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human Tetraspan 5
Isotype:	IgG
Predicted Reactivity:	Human, Mouse, Rat, Dog, Cow, Pig, Horse, Chicken, Rabbit
Purification:	Purified by Protein A.

### Target Details

Target:	Tetraspanin 5 (TSPAN5)
Alternative Name:	Tetraspan 5 ( <a href="#">TSPAN5 Products</a> )
Background:	Synonyms: NET 4, TSPAN 5, NET4, Tetraspan 5, Tetraspan-5, Tetraspan NET 4, Tetraspan

## Target Details

---

TM4SF, Tetraspanin 5, Tetraspanin5, Tetraspan5, Tspaw5, TM4SF9, Transmembrane 4 superfamily member 8, Transmembrane 4 superfamily member 9.

Background: TSPAN5 is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility.

---

Gene ID: 10098

## Application Details

---

Application Notes: IF(IHC-P) 1:50-200  
IF(IHC-F) 1:50-200  
IF(ICC) 1:50-200

---

Restrictions: For Research Use only

## Handling

---

Format: Liquid

---

Concentration: 1 µg/µL

---

Buffer: Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.

---

Preservative: ProClin

---

Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.

---

Storage: -20 °C

---

Storage Comment: Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.

---

Expiry Date: 12 months