

Datasheet for ABIN708693

anti-APOA2 antibody (AA 24-100) (FITC)[Go to Product page](#)

Overview

Quantity:	100 µL
Target:	APOA2
Binding Specificity:	AA 24-100
Reactivity:	Human, Rat
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This APOA2 antibody is conjugated to FITC
Application:	Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human APOA2
Isotype:	IgG
Cross-Reactivity:	Human, Rat
Predicted Reactivity:	Mouse,Dog
Purification:	Purified by Protein A.

Target Details

Target:	APOA2
Alternative Name:	APOA2 (APOA2 Products)

Target Details

Background:	Synonyms: apoAII, Apo-AII, ApoA-II, Apolipoprotein A-II, Apolipoprotein A2, APOA2 Background: Apolipoprotein A II is the second most abundant protein of the high density lipoprotein particles. The apolipoprotein A II gene consists of 4 exons and 3 introns. The four exons encode the 5' untranslated region, pre peptide, a short N terminal domain and a C terminal domain composed of a variable number of lipid binding amphipathic helices. Familial apolipoprotein A II deficiency may result from a splice junction alteration which blocks splicing of intron 3 from the primary transcript and result in the formation of a non functional mRNA.
Gene ID:	336
UniProt:	P02652
Pathways:	Regulation of Lipid Metabolism by PPARalpha , Production of Molecular Mediator of Immune Response , Negative Regulation of Transporter Activity , Lipid Metabolism

Application Details

Application Notes:	IF(IHC-P) 1:50-200 IF(IHC-F) 1:50-200 IF(ICC) 1:50-200
Restrictions:	For Research Use only

Handling

Format:	Liquid
Concentration:	1 µg/µL
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	-20 °C
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.
Expiry Date:	12 months