antibodies

## Datasheet for ABIN7091234 anti-HPV16 E7 antibody (Alexa Fluor 680)



Overview

Quantity:	100 μL
Target:	HPV16 E7 (HPV-16 E7)
Reactivity:	Human Papillomavirus 16 (HPV-16), Virus
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This HPV16 E7 antibody is conjugated to Alexa Fluor 680
Application:	Western Blotting (WB), Immunofluorescence (Paraffin-embedded Sections) (IF (p))
Product Details	
Immunogen:	KLH conjugated synthetic peptide derived from HPV16 E7 protein
lsotype:	lgG
Cross-Reactivity:	Virus
Cross-Reactivity (Details):	HPV16
Purification:	Purified by Protein A.
Target Details	

Target:	HPV16 E7 (HPV-16 E7)
Alternative Name:	HPV16 E7 (HPV-16 E7 Products)
Target Type:	Viral Protein
Background:	Synonyms: E7, Human Papilloma Virus, Human papillomavirus type 16 E7, Human

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN7091234 | 03/08/2024 | Copyright antibodies-online. All rights reserved.

	papillomavirus type 16 E7, Protein E7, Human Papillomavirus 16 E7, HPV16 E7, E7 protein
	HPV16, VE7_HPV16.
	Background: Human papilloma viruses (HPVs) can be classified as either high risk or low risk
	according to their association with cancer. HPV16 and HPV18 are the most common of the
	high risk group while HPV6 and HPV11 are among the low risk types. Approximately 90 $\%$ of
	cervical cancers contain HPV DNA of the high risk types. Mutational analysis have shown that
	the E6 and E7 genes of the high risk HPVs are necessary and sufficient for HPV transforming
	function. The specific interactions of the E6 and E7 proteins with p53 and pRB, respectively,
	correlate with HPV high and low risk classifications. The high risk HPV E7 proteins bind to pRB
	with a higher affinity than do the low risk HPV proteins, and only the high risk HPV E6 proteins
	form detectable complexes with p53 in vitro.
Gene ID:	1489079
Application Details	
Application Notes:	IF(IHC-P) 1:50-200
Restrictions:	For Research Use only
Handling	
Handling Concentration:	1 μg/μL
	1 μg/μL Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and
Concentration:	
Concentration:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and
Concentration: Buffer:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Concentration: Buffer: Preservative:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol. ProClin
Concentration: Buffer: Preservative:	Aqueous buffered solution containing 0.01M TBS ( pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.   ProClin   This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be