antibodies -online.com







anti-BBS7 antibody



| () | ve | K\ / | | A . |
|-----|--------|------|----|-----|
| | \cup | 1 V/ | Щ. | V۷ |
| | | | | |

| Quantity: | 100 μg |
|--------------|--|
| Target: | BBS7 |
| Reactivity: | Human, Mouse |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Conjugate: | This BBS7 antibody is un-conjugated |
| Application: | ELISA, Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP) |

Product Details

| Immunogen: | Bardet-Biedl syndrome 7 |
|---------------|---------------------------------|
| Isotype: | IgG |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95 % as determined by SDS-PAGE |

Target Details

| Target: | BBS7 | |
|-------------------|--|--|
| Alternative Name: | BBS7 (BBS7 Products) | |
| Background: | Synonyms:BBS2L1 Background:The BBSome complex is thought to function as a coat complex required for sorting of specific membrane proteins to the primary cilia. The BBSome complex is required for ciliogenesis but is dispensable for centriolar satellite function. This ciliogenic function is mediated in part by the Rab8 GDP/GTP exchange factor, which localizes to the basal | |

Target Details

body and contacts the BBSome. Rab8(GTP) enters the primary cilium and promotes extension of the ciliary membrane. Firstly the BBSome associates with the ciliary membrane and binds to RAB3IP/Rabin8, the guanosyl exchange factor(GEF) for Rab8 and then the Rab8-GTP localizes to the cilium and promotes docking and fusion of carrier vesicles to the base of the ciliary membrane. The BBSome complex, together with the LTZL1, controls SMO ciliary trafficking and contributes to the sonic hedgehog(SHH) pathway regulation. Required for proper BBSome complex assembly and its ciliary localization.

Molecular Weight: 80 kDa

Gene ID: 55212

UniProt: Q8IWZ6

Pathways: Hedgehog Signaling

Application Details

Application Notes: WB: 1:500-1:3000, IP: 1:500-1:1000, IHC: 1:20-1:200

Restrictions: For Research Use only

Handling

| Format: | Liquid | |
|--------------------|--|--|
| Buffer: | PBS with 0.02 % sodium azide and 50 % glycerol pH 7.3, | |
| Preservative: | Sodium azide | |
| Precaution of Use: | This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only. | |
| Storage: | -20 °C | |
| Storage Comment: | -20°C for 12 months (Avoid repeated freeze / thaw cycles.) | |
| Expiry Date: | 12 months | |