



[Go to Product page](#)

Datasheet for ABIN7154199
anti-PAEP antibody (AA 19-180) (Biotin)

Overview

| | |
|----------------------|--|
| Quantity: | 100 µg |
| Target: | PAEP |
| Binding Specificity: | AA 19-180 |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Conjugate: | This PAEP antibody is conjugated to Biotin |
| Application: | ELISA |

Product Details

| | |
|-------------------|---|
| Immunogen: | Recombinant Human Glycodelin protein (19-180AA) |
| Isotype: | IgG |
| Cross-Reactivity: | Human |
| Purification: | >95%, Protein G purified |

Target Details

| | |
|-------------------|---|
| Target: | PAEP |
| Alternative Name: | PAEP (PAEP Products) |
| Background: | Background: This protein is, quantitatively, the main protein synthesized and secreted in the endometrium from mid-luteal phase of the menstrual cycle and during the first semester of |

Target Details

pregnancy.

Aliases: Alpha uterine protein antibody, gD antibody, GdA antibody, GdF antibody, GdS antibody, Glycodelin A antibody, Glycodelin antibody, Glycodelin F antibody, Glycodelin S antibody, MGC138509 antibody, MGC142288 antibody, PAEG antibody, PAEP antibody, PAEP_HUMAN antibody, PEG antibody, PEP antibody, Placental Protein 14 / Glycodelin A antibody, Placental protein 14 antibody, PP14 antibody, Pregnancy-associated endometrial alpha-2 globulin antibody, Progesterone-associated endometrial protein antibody, Progesterone-associated endometrial protein antibody

UniProt: [P09466](#)

Application Details

Application Notes: Optimal working dilution should be determined by the investigator.

Restrictions: For Research Use only

Handling

Format: Liquid

Buffer: Preservative: 0.03 % Proclin 300
Constituents: 50 % Glycerol, 0.01M PBS, pH 7.4

Preservative: ProClin

Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.

Storage: -20 °C,-80 °C

Storage Comment: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.