

Datasheet for ABIN719225 anti-El24 antibody (AA 251-340) (HRP)



Overview

Background:

Overview	
Quantity:	100 μL
Target:	El24
Binding Specificity:	AA 251-340
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This EI24 antibody is conjugated to HRP
Application:	ELISA, Immunohistochemistry (Paraffin-embedded Sections) (IHC (p)), Immunohistochemistry (Frozen Sections) (IHC (fro))
Product Details	
Immunogen:	KLH conjugated synthetic peptide derived from human TP53I8/EI24
Isotype:	IgG
Predicted Reactivity:	Human,Mouse,Rat,Dog,Cow,Pig,Horse
Purification:	Purified by Protein A.
Target Details	
Target:	El24
Alternative Name:	EI24 (EI24 Products)
Da alamana da	C

Synonyms: Ei24, EI24_HUMAN, EPG4, Etoposide induced 2.4 mRNA, Etoposide-induced protein

2.4 homolog, p53-induced gene 8 protein, PIG8, TP53I8, Tumor protein p53 inducible protein 8. Background: This gene has higher expression in p53-expressing cells than in control cells and is an immediate-early induction target of p53-mediated apoptosis. The protein encoded by this gene contains six putative transmembrane domains and may suppress cell growth by inducing apoptotic cell death through the caspase 9 and mitochondrial pathways. This gene is located on human chromosome 11q24, a region frequently altered in cancers. Alternative splicing results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008].

Gene ID:

9538

Application Details

Restrictions:	For Research Use only
	IHC-F 1:100-500
Application Notes:	IHC-P 1:200-400

Handling

Format:	Liquid
Concentration:	1 μg/μL
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Handling Advice:	Do NOT add Sodium Azide! Use of Sodium Azide will inhibit enzyme activity of horseradish peroxidase.
Storage:	-20 °C
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.
Expiry Date:	12 months