

Datasheet for ABIN7427114 **anti-Ephrin B2 antibody**



[Go to Product page](#)

Overview

Quantity:	100 µL
Target:	Ephrin B2 (EFNB2)
Reactivity:	Human
Host:	Mouse
Clonality:	Monoclonal
Conjugate:	This Ephrin B2 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP), Immunocytochemistry (ICC)

Product Details

Purpose:	Monoclonal Antibody to Ephrin B2 (EFNB2)
Immunogen:	Recombinant Ephrin B2 (EFNB2) corresponding to N-terminal His and GST Tag
Isotype:	IgG
Specificity:	The antibody is a mouse monoclonal antibody raised against EFNB2. It has been selected for its ability to recognize EFNB2 in immunohistochemical staining and western blotting.
Purification:	Protein A + Protein G affinity chromatography

Target Details

Target:	Ephrin B2 (EFNB2)
Alternative Name:	EFNB2 (EFNB2 Products)

Target Details

Background:	EPLG5, HTKL, Htk-L, LERK5, HTK Ligand, Ligand Of Eph-Related Kinase 5, Eph-Related Receptor Tyrosine Kinase Ligand 5
UniProt:	P52799
Pathways:	RTK Signaling, Regulation of Muscle Cell Differentiation

Application Details

Application Notes:	Western blotting: 0.5-2 µg/mL, Immunohistochemistry: 5-20 µg/mL, Immunocytochemistry: 5-20 µg/mL, Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.
Restrictions:	For Research Use only

Handling

Format:	Liquid
Concentration:	1 mg/mL
Buffer:	PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.
Preservative:	Sodium azide
Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C, -20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Expiry Date:	24 months