

Datasheet for ABIN7431143
anti-TIMP4 antibody (AA 50-224)[Go to Product page](#)

1 Image

Overview

Quantity:	100 µL
Target:	TIMP4
Binding Specificity:	AA 50-224
Reactivity:	Cow
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This TIMP4 antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunoprecipitation (IP), Immunocytochemistry (ICC)

Product Details

Purpose:	Polyclonal Antibody to Tissue Inhibitors Of Metalloproteinase 4 (TIMP4)
Immunogen:	Recombinant Tissue Inhibitors Of Metalloproteinase 4 (TIMP4) corresponding to Ala50~Pro224 (Accession # 097563)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against TIMP4. It has been selected for its ability to recognize TIMP4 in immunohistochemical staining and western blotting.
Cross-Reactivity:	Human, Mouse, Rat
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography

Target Details


Target:	TIMP4
Alternative Name:	Tissue Inhibitors Of Metalloproteinase 4 (TIMP4 Products)
Background:	TIMP Metallopeptidase Inhibitor 4, Tissue Inhibitor of Metalloproteinases 4

Application Details

Application Notes:	Western blotting: 0.5-2 µg/mL 1:200-800 Immunohistochemistry: 5-20 µg/mL 1:20-80 Immunocytochemistry: 5-20 µg/mL 1:20-80 Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.
Restrictions:	For Research Use only

Handling

Format:	Liquid
Buffer:	PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.
Preservative:	Sodium azide
Precaution of Use:	This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Expiry Date:	24 months



Western Blotting

Image 1. Detection of Recombinant TIMP4, Bovine using Polyclonal Antibody to Tissue Inhibitors Of Metalloproteinase 4 (TIMP4)