



[Go to Product page](#)

Datasheet for ABIN7439825
anti-S-Adenosyl-Methionine antibody

Overview

Quantity:	100 µL
Target:	S-Adenosyl-Methionine (SAM)
Reactivity:	Various Species
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This S-Adenosyl-Methionine antibody is un-conjugated
Application:	Chemiluminescence Immunoassay (CLIA), ELISA, Immunocytochemistry (ICC), Immunohistochemistry (Frozen Sections) (IHC (fro))

Product Details

Purpose:	Polyclonal Antibody to S-Adenosyl Methionine (SAM)
Immunogen:	OVA Conjugated S-Adenosyl Methionine (SAM)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against SAM. It has been selected for its ability to recognize SAM in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography followed by Protein A affinity chromatography

Target Details

Target:	S-Adenosyl-Methionine (SAM)
Alternative Name:	S-Adenosyl Methionine (SAM Products)

Target Details

Target Type:	Chemical
Background:	SAMe, SAM-e, S-Adenosylmethionine, S-(5'-Adenosyl)-L-Methionine Chloride

Application Details

Application Notes:	Enzyme-Linked Immune Sorbent Assay: 212.5 ng/mL Immunohistochemistry: 5-20 µg/mL Immunocytochemistry: 5-20 µg/mL Optimal working dilutions must be determined by end user.
Comment:	The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.
Restrictions:	For Research Use only

Handling

Format:	Liquid
Buffer:	0.01M PBS, pH 7.4, containing 0.05 % Proclin-300, 50 % glycerol.
Preservative:	ProClin
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.
Storage:	4 °C,-20 °C
Storage Comment:	Store at 4°C for frequent use. Stored at -20°C in a manual defrost freezer for two year without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Expiry Date:	24 months